# AbraPlan-10

## Instruction Manual

Manual No.: 15497001

Date of Release 16.03.2005

*AbraPlan-10*  
*Instruction Manual*

*AbraPlan-10*  
*Instruction Manual*

<b>Table of Contents</b>	<b>Page</b>
User's Guide .....	1
Reference Guide .....	25
Quick Reference Guide.....	36

---

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

**Instruction Manuals:** Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

**Service Manuals:** Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2005.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup  
Denmark  
Telephone +45 44 600 800  
Fax +45 44 600 801

---

## **AbraPlan-10**

### **Safety Precaution Sheet**

#### **To be read carefully before use**

1. The operator should be fully aware of the use of the machine according to the Instruction Manual.
2. The machine must be installed in compliance with local safety regulations.
3. Be aware that the machine's centre of gravity is located in the upper half of the machine.
4. Before lifting the machine, ensure that the supplied lifting boom is securely bolted to the machine.
5. When lifting the machine using a forklift, lift from front or rear - never lift the machine from the side.
6. When lifting the machine using lifting straps, ensure that the straps are crossed and do not press on the sides of the machine.
7. The machine must be placed on a safe and stable support, which is capable of bearing the weight of this machine. Before using the machine, it must be levelled using the adjustable legs.
8. Be sure that the actual voltage corresponds to the voltage stated on the side of the machine and that the connections comply with local regulations. The machine must be earthed.
9. Be aware that when the machine is connected to a compressed air supply the specimen holder arm moves upwards.
10. Make sure that the specimens in the specimen holder are securely fixed, and ensure that the securing screws are not sticking out.
11. Keep clear of the grinding stone when operating.

*AbraPlan-10*  
*Instruction Manual*

- 12.** Always use Safety Goggles when operating the machine.
- 13.** If you observe malfunctions or hear unusual noises - stop the machine and call technical service.
- 14.** The machine must be disconnected from the mains supply prior to any service and before connecting or disconnecting the recirculation pump.
- 15.** To achieve maximum safety and lifetime of the machine, use only original Struers consumables.

---

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

---

# User's Guide

<b>Table of Contents</b>	<b>Page</b>
<b>1. Getting Started</b>	
Checking the Contents of the Packing Box .....	2
AbraPlan-10.....	2
Recirculation Cooling Unit (optional).....	2
Placing AbraPlan-10.....	2
Getting Acquainted with AbraPlan-10 .....	3
Noise Level.....	3
Supplying Power .....	4
Direction of Rotation .....	4
Supplying Compressed Air.....	4
Connecting the Recirculation Cooling Unit (accessory).....	5
Mounting the Grinding Stone/ Diamond Grinding Disc.....	6
Flushing Unit Head.....	8
Mounting Stock Removal Sensor (accessory) .....	9
Mounting the Water Level Indicator (accessory) .....	9
<b>2. Basic Operations</b>	
Front Panel .....	10
Front Panel Controls .....	10
Software Settings.....	11
Setting the Language .....	12
Reading the Display .....	14
Sleep Mode.....	15
Changing/Editing Values .....	15
Numeric Values .....	15
Text Values .....	16
Setting up the Software.....	17
Dressing Arm Movement .....	18
Inserting/Removing the Specimen Holder .....	18
Inserting the Specimen Holder.....	18
Removing the Specimen Holder.....	19
Grinding Setup .....	19
Setting the Process Time .....	20
Setting Stock Removal (accessory) .....	20
Setting the Force .....	21
Cooling Water .....	21
Starting the Preparation Process .....	22
Stopping the Preparation Process.....	22
Dressing the Grinding Stone.....	23
Adjusting the Position of the Dressing Arm .....	23
Dressing the Diamond Grinding Disc.....	24

## **1. Getting Started**

### **Checking the Contents of the Packing Box**

*AbraPlan-10*

In the packing box you should find the following parts:

- 1 AbraPlan-10
- 1 Outlet hose (factory mounted)
- 1 Inlet hose (factory mounted)
- 2 Hose clamps, 17 mm
- 1 Drain angle, 90°
- 1 Hose clamp, 35-60 mm
- 1 Hose for compressed air
- 2 Hose clamps, 12 mm
- 1 Hose connection for compressed air
- 1 Rubber disc 350 mm
- 1 Flange
- 1 Bolt M12
- 1 Allen key 8 mm
- 1 Set of Instruction Manuals

*Recirculation Cooling Unit*  
*(optional)*

- 1 Recirculation Cooling Unit
- 1 Drain angle, 45°
- 1 Funnel
- 2 Plastic inserts

### **Placing AbraPlan-10**


AbraPlan-10 should be placed directly on a plane and horizontal floor, capable of bearing the weight of this machine, please refer to the Technical Data section.

**Getting Acquainted with  
 AbraPlan-10**


Take a moment to familiarise yourself with the location and names of the AbraPlan-10 components.

**Front view**

- ① Shield
- ② Control panel,  
(see 2. Basic Operations)
- ③ Specimen holder motor
- ④ Specimen quick coupling
- ⑤ Emergency Stop
- ⑥ Cover plate for  
Recirculation Cooling Unit

**Side views**


- ⑦ Compressed air inlet
- ⑧ Cast iron beam
- ⑨ Mains power switch
- ⑩ Two-handed start buttons
- ⑪ Bottom cover
- ⑫ Electrical connections

**Noise Level**

Approx. 68 dB (A) measured during idle running, at the operator's position in front of the machine.


## Supplying Power

### **IMPORTANT**

Check that the mains supply voltage corresponds to the voltage stated on the Type Plate (located under the mains switch on the side of the machine).  
If the machine is already connected to a mains supply, disconnect this supply before removing the Allen screws securing the electrical panel.

- Remove the 4 Allen screws, on the right hand side of the machine, securing the electrical panel, and let the electrical panel rest on the tabs.
- Lead the cable through the conduit in the panel and connect the 3 phases and earth according to local regulations.
- Where an external recirculation unit is used, connect the electrical cable from the pump, according to the diagram located inside the panel.

## Direction of Rotation

Check that when the power is turned on the grinding stone/ diamond grinding disc rotates counter-clockwise.  
If this is not the case, switch off AbraPlan-10 and unplug the machine.


- Change two of the phases.
- Repeat the rotation check.

## Supplying Compressed Air

- Connect the compressed air supply to the inlet, located at the rear of the left-hand side of the machine, using the air hose and the hose connection delivered with the machine.
- Fasten the air hose with a hose clamp.

The pressure supply should be 6-10 bar and can be supplied either from a central compressor, a portable compressor with a compressed air reservoir, or a compressed-air bottle. A capacity of 20 l/min at atmospheric pressure is sufficient.

## Connecting the Recirculation Cooling Unit (accessory)

### **Warning**

Disconnect the machine from the mains power supply before connecting or disconnecting the recirculation pump.

- Place the unit where you find it convenient, either inside the machine or outside, on the left.
- Slide the metal spring over the inlet hose (this will prevent kinks in the hose that could obstruct the flow of water).
- Connect the inlet hose from the flushing unit to the pump of the recirculation unit.
- Secure the metal spring in place by turning two coils around the connector (see illustration).
- Cut the drain hose to the correct length, remove approximately 20 mm of the spring inside the drain hose and fit the hose over the 90° drain angle.
- Slide the drain angle over the pipe in the top of the reservoir lid. (Note: The 45° elbow pipe delivered with the Recirculation Cooling Unit is not used with this machine).
- Check that there is a steady fall on the whole course of the outlet hose.
- Open the electrical connections cover
  - If mounting inside the machine, lead the cable into the electrical panel from inside the machine through the rear.
  - If mounting outside the machine, lead the cable through the conduit in the panel.
- Now connect the cable to the AbraPlan-10 using the terminal connectors in the electrical panel.
- Check that when the power is turned on, the pump rotates in the direction indicated by the arrow on the top of the pump. If this is not the case, change 2 of the electrical phases.
- Replace the cover to the electrical connections.
- If mounting inside the machine, remove the front cover and place the tank on the trolley.
- Place a disposable plastic insert in the tank and fold it over the rim.
- Fill the tank with 29.1 l water and 900 ml Struers Additive for Cooling Fluid.  
The water level should be 8-10 cm below the upper rim of the tank.

### **IMPORTANT**

Too high a level of coolant in the tank might damage the pump. To avoid this, place the disposable insert so that the pre-punched hole is in front of the overflow aperture in the tank.

**IMPORTANT**


Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the Additive container). Remember to top up with Struers Additive each time you refill with water.

- Fit the tank lid and sieve.
- Move the unit into place and lead the outlet and overflow drain hoses into the cooling water tank. Adjust the lengths if necessary.
- If mounting inside, replace the front cover.

**Mounting the Grinding Stone/  
Diamond Grinding Disc**

**IMPORTANT**


Make sure that the grinding stone/diamond grinding disc is intact. The stone/disc must be dry when mounted and the flange should be clean and smooth.

- ① Flushing unit head
- ② Securing finger screws
- ③ Stone guard
- ④ Dressing arm

- Ensure that the sample motor is fully raised, and lift the shield to gain access to the grinding area.
- Make sure that the dressing arm is located in its storage position (on the left-hand side).
- Lift the flushing unit head and, using its built-in magnet, attach it to the sample motor mounting arm.
- Move the sliding cover on the stone guard across until it is free at both ends.
- Unscrew the 2 finger screws (located on the right-hand side)
- Use both the handle and the flushing unit connection to lift the stone guard up and away to the right.

- Assemble the grinding stone/ diamond grinding disc on the flange as illustrated:

- ① Motor Flange
- ② Rubber disk
- ③ Grinding stone/ diamond grinding disc
- ④ Securing flange and cardboard washer
- ⑤ Securing bolt

***IMPORTANT***


Do not over tighten the securing bolt as this may damage the grinding stone/disc.

- Mount the bolt and fasten firmly.
- Replace the stone guard and secure it with the 2 finger screws.

***IMPORTANT***

Remember to re-attach the flushing unit head.  
Remember to raise the dressing arm to the top of its travel.

## Flushing Unit Head

When attached to its mounting (as illustrated), the flushing unit supplies water/coolant to the grinding stone/ diamond grinding disc during the grinding and dressing processes. However, the flushing unit can be removed from the mounting and used to hose down the grinding area and samples etc. To do this:

- Squeeze in and hold the clamp buttons (illustrated by arrows) to cut-off the water/coolant flow.
- Press F1 to start the pump.
- Lift the flushing unit free of its mounting.
- Direct the flushing unit in the desired direction and release the clamp buttons.
- Press F1 to stop the pump.
- Reinsert the flushing unit in its mounting.


If the grinding/dressing process is not active, a water/coolant flow can be achieved by pressing F1 when the GRINDING or DRESSING Menus are displayed.

### Mounting the Stock Removal Sensor (accessory)


- Disconnect main power.
- Connect the plug on the removal sensor with the plug located on the inside of the control panel.
- Take the removal sensor and gently slide it up in the hole located in the bottom of the cast iron beam.
- Let the bracket slide over the factory installed bracket.
- From the left, slide the two screws through the brackets, and tighten them with a 10 mm spanner.
- Reconnect the main power and follow the instructions on the display, using the tools supplied and the calibrating sticks.

### Mounting the Water Level Sensor (accessory)


- Switch the machine off.
- Pull out the recirculation unit.
- Remove the pump, drain hose and the top cover.
- Remove the red cover from the hole for the sensor or remove the existing Water Level Indicator.
- Remove the locking nut on the bottom of the sensor and slide off the black float.
- Mount the new Water Level Sensor on the recirculation unit cover and secure using the nut supplied (flat side of nut up).
- Replace the float (with the clear surface up) and secure the locking nut in place.
- On the right hand side of the machine, remove the 4 Allen screws securing the electrical panel, and let the electrical panel hang on the two tabs (as illustrated in Supplying Power).
- Remove the dummy jack plug located on the electrical panel by pressing down the red tab.
- Plug in the jack plug from the sensor.
- Reassemble the recirculation unit and slide it back into position.
- Replace the front cover.
- Switch the machine on.






## 2. Basic Operations

### Front Panel


### Front Panel Controls

Name	Key	Function	Name	Key	Function
FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	Push/Turn Knob		Used for entering and changing parameters. Combined cursor and enter key.
FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	Esc		Leaves the present menu or aborts functions/changes.
FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	STOP		Stops the preparation process.
FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	EMERGENCY STOP		The EMERGENCY STOP is located on the front of the machine. - Push the red button to stop. - Pull the red button to release.
			MAIN SWITCH		The main switch is located on the right side of the machine.

Switch on the power at the main switch located at the right hand side of the machine. The following display will appear briefly:

## Software Settings

The display will then change to the last screen shown before AbraPlan-10 was switched off. When switching on AbraPlan-10 for the first time, the MAIN MENU display will appear. If the heading in the display is different, press **Esc**, until the MAIN MENU appears.

The MAIN MENU is the highest level in the menu structure. From this menu, you can enter the configuration menu and grinding process menu.

*Setting the Language*


Before continuing, we recommend that you select the language that best suits you.

Turn knob to select CONFIGURATION

Push knob to activate the CONFIGURATION Menu

Turn knob to select Language.

Push knob to activate the Language pop-up menu.

Turn knob to select the language you prefer.

*AbraPlan-10*  
*Instruction Manual*

Push knob to accept the language.


The CONFIGURATION Menu now appears in the language you have chosen.

**Esc**

Press **Esc** to return to the MAIN MENU.


**Reading the Display**

The display is primarily divided into 4 areas. The position of these areas and the information they contain are explained in the illustration below, which uses the CONFIGURATION Menu as an example:

- A** Heading: this tells you where you are in the software.
- B** Information fields: these will either be numerical values or text fields providing information associated with the process shown in the heading. The inverted text shows the cursor position.
- C** Help text: helpful hints to allow you to continue.
- D** Function key options: the functions of these change with the window displayed.

During the grinding process the screen could look as follows:

- A** Grinding time remaining
- B** Force applied on specimen holder
- C** Load on main motor
- D** Function(s) selectable during process

*Sleep Mode*

To increase the lifetime of the display, the backlight is automatically switched off if AbraPlan-10 is not used for 15 min. Press any key to re-activate the backlight.

**Please Note**

The sample screens in this Instruction Manual show a number of possible texts. The actual screen displayed may differ from the samples in this manual.


**Changing/Editing Values**

*Numeric Values*

Depending on the type of value, there are two different ways of editing.


Turn knob to select the value to be changed, e.g. Process time:

Push knob to edit value.

Two square brackets [ ] appear around the value.

Turn knob clockwise to increase, or counter-clockwise to decrease the numeric value.


Push knob to accept the new value. (Pressing **Esc**, aborts the changes, preserving the original value.)

*Text Values*

Turn knob to select the text value to be changed, e.g. Language

Push knob to edit the value.

If there are only two options, then a toggle function is active, e.g. Yes / No.

If there are more than two options, a pop-up menu appears:

Turn knob to select/toggle the correct option.

Push knob to accept the new value. (Pressing **Esc**, aborts the changes, preserving the original setting.)

## Setting up the Software


Before operating the AbraPlan-10, you should check/edit the software configuration values:

Turn knob to select CONFIGURATION.

Push knob to activate the CONFIGURATION Menu.

Turn knob to select desired parameter.

Push knob to edit the parameter setting.

The settings possible and the default settings are shown below:

	Possibilities	Increment	Default
Display contrast	0 - 50	1	25
Language	English German French Spanish Japanese		English
Water when dressing	Yes/No		Yes
Dressing arm movement (mm)	0.02 – 0.50	0.02	0.10

Repeat the process until all values are correct.

Press **Esc** to return to the MAIN MENU.


*Dressing Arm Movement*

This setting determines the distance the dressing arm is lowered each time the dressing arm handle is lifted or the raise and lower function keys are pressed when in the DRESSING menu.

**IMPORTANT**

Do not repeatedly lift the dressing arm handle between grinding stone/diamond grinding disc dressings as this may result in excessive removal of the stone/disc when dressing the next time.

**Inserting/Removing  
the Specimen Holder**

*Inserting the Specimen Holder*

- Position the specimen holder under the quick coupling, and support it with your fingertips.
- Press and hold the flange of the column down with the heel of your hand while guiding the pressure tap of the specimen holder into the coupling.
- Turn the specimen holder until the three pins engage with the corresponding holes. Check that the specimen holder is in position by trying to turn it to the side.
- Release the flange.


*Removing the Specimen Holder*

- Place your hands as shown in the illustration above.
- Press upwards with your fingers to lift the specimen holder slightly. At the same time, press and hold down the flange with the heel of your hand.
- Using your fingers to support the specimen holder; lower it free of the coupling.
- Release the flange and completely remove the specimen holder.

**Grinding Setup**

Before starting the grinding process, the grinding time and force setting should be checked/changed. This is done using the control panel software:

- ↓ If the MAIN MENU is not currently displayed, press **Esc** twice
-  Turn knob to select GRINDING:


- ↓
-  Push knob to enter the GRINDING MENU.

↓

*Setting the Process Time*

Turn knob to select Process time setting.

Push knob to display square brackets [ ] and turn knob to change the Process time setting.

Note that stock removal is not available when process time is selected.


Push knob to accept the new value.

*Setting Stock Removal  
(accessory)*

Turn knob to select Stock removal setting.

Push knob to display square brackets [ ] and turn knob to change the Stock removal setting.

↓ Note: the Process time is automatically disabled.

Push knob to accept the new value.


The Stock removal process has a 15 minutes time-out. If the prescribed amount of material has not been removed within this time, the process stops automatically.

### *Setting the Force*

Turn knob to select Force setting.

Push knob to display square brackets [ ] and turn knob to change the Force setting.

Push knob to accept the new value.

### *Cooling Water*

The cooling water pump will automatically start when the preparation process is started. To stop the pump, press F1.


**Please Note**

All the value and function settings defined in the GRINDING SETUP menu are saved in the AbraPlan-10's memory. Therefore, if the power to the machine is interrupted, these values are remembered.

**Starting the  
Preparation Process**

- Insert the specimen holder.
- Ensure that the dressing arm is located in the “parked or resting” position on the left-hand side of the machine and that the sliding cover on the stone guard is located all the way to the right.
- Lower the shield.
- Using the software display, enter the GRINDING menu and set/check the correct time and force.
- Start the grinding process by pressing and holding the two start buttons simultaneously until the specimen holder starts rotating.

The software display will appear as:

**Please Note**

Although the force setting can be changed during the grinding process, this is not recommended when using Stock removal.

**Stopping the  
Preparation Process**

When the time has elapsed, the grinding stone will automatically stop rotating and the specimen holder will return to its initial position.

## Dressing the Grinding Stone

### **IMPORTANT**

The grinding process can be stopped at any time by pressing the STOP button on the Control Panel or by pressing the EMERGENCY STOP.

AbraPlan-10 is fitted with a diamond tool for dressing the grinding stone. It is important to dress the grinding stone at regular intervals to keep the stone plane and sharp.

### **Important**

The dressing arm with the diamond tool should never be used on a diamond grinding stone since this would destroy the disc or the diamond tool.


When a new grinding stone is fitted, the AbraPlan-10 software asks by default if you want to raise the dressing arm to the top of its movement.

The dressing arm moves down automatically each time the handle on the dressing arm is lifted. The distance of this incremental arm movement is defined by the “Dress. arm movement:” value in the CONFIGURATION menu.

The position of the dressing arm can be checked and changed through the software display on the Control Panel via the DRESSING menu.

### *Adjusting the Position of the Dressing Arm*

When the sliding cover on the stone guard is covering the grinding stone (i.e. moved all the way over to the left) and the dressing arm handle is lifted, the software automatically displays the DRESSING menu.

The DRESSING menu shows the current height or “Position” of the dressing arm. This position can be changed using F3 to raise the dressing arm and F4 to lower it.

Each time F3 or F4 is pressed, the dressing arm will move in the direction specified by the incremental distance defined in the CONFIGURATION menu.

To move the dressing arm more than the defined incremental distance, press the function repeatedly or press and hold to move the arm to end of its travel (either top or bottom).

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↓ Check the position of the dressing arm.

**F3** If the position is incorrect, press F3 to raise the dressing arm.

Or

**F4** Press F4 to lower the dressing arm.

↓

- Swing the diamond tool to a position over the grinding stone and move the diamond tool 2-6 times forwards and backwards across the grinding stone in a slow and even movement without pressing against the stone. If necessary, adjust the position of the diamond tool as described above.
- When the dressing is complete, swing the diamond tool back to its “parked or resting” position.

The software display returns to the menu displayed immediately prior to the dressing.

## **Dressing the Diamond Grinding Disc**

In order to dress a diamond grinding stone, mount 3 aluminium oxide dressing sticks in a sample holder and grind for a few seconds.

# Reference Guide

Table of Contents	Page
<b>1. Accessories and Consumables</b> .....	26
<b>2. Struers Metalog Guide™</b> .....	27
<b>3. Trouble-Shooting</b> .....	28
<b>4. Maintenance</b>	
Daily Service .....	32
Weekly Service .....	32
Checking the Recirculation Cooling Unit.....	32
Monthly Service .....	32
Replacing the Cooling Water .....	32
Adjusting the V-belts .....	33
<b>5. Technical Data</b> .....	34

## 1. Accessories and Consumables

Specification	Code
<i>Grinding stone</i> , Al <sub>2</sub> O <sub>3</sub> , for plane grinding of soft and ductile metals (<HV 250), 356 mm dia., grit 60	ABGOF
<i>Grinding stone</i> , Al <sub>2</sub> O <sub>3</sub> , for plane grinding of harder metals (>HV 250), 356 mm dia., grit 150	ABGAL
<i>Grinding stone</i> , SiC, for plane grinding of non-ferrous metals, 356 mm dia., grit 150	ABGAS
<i>Diamond grinding disc</i> , 356 mm dia., grit 120	ABWEE
<i>Diamond point tool</i>	PLAMT
<i>Accurate Material Removal Unit</i> Stops process after required amount of material has been removed ( <i>Stock Removal Sensor</i> )	ABPEU
<i>Cooling unit</i> Recirculation Cooling Unit, complete, 3-phased	TRECI
<i>Rolling Pallet for TRECI</i>	TREJA
<i>Water Level Indicator</i> Sensor indicates the level of water in the Recirculation Cooling Unit (TRECI)	PAMWA
<i>Additive for Cooling Fluid</i> 1 l 5 l	ADDUN ADDFI
<i>Bacterial Cleaning Liquid</i> 0.5 l	BACTE
<i>Disposable Plastic Insert</i> for collecting waste material in the Recirculation Cooling Unit. Package with 20 pcs.	TREPO

***Remember...***

Struers offers a comprehensive range of consumables for grinding and polishing. Please ask for separate leaflets.

## **2. Struers Metalog Guide™**

In Struers Metalog Guide™ you will find a detailed description of grinding/polishing methods for automated mechanical specimen preparation.

Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method is easy, including choice of consumables. Always consult Struers Metalog Guide™ for the correct preparation method for the actual specimens.

Struers Metalog Guide™ contains 6 useful chapters:

Metalogram: a quick and safe guide to the right preparation method.

Metalog Methods: a complete catalogue of preparation methods, based on Struers' vast store of materialographic experience, and employing Struers' range of consumables.

Preparation Philosophy: the basics of modern specimen preparation, seen from a professional point of view.

Metalog Process: the materialographic preparation process from start to finish, logically explained.

Metalog Master: a combined trouble-shooting guide and supply of in-depth information on the processes of mechanical preparation, including an expert system for the solving of preparation problems.

Consumables Specification: quick access to the relevant consumables for the chosen preparation methods.

A complete guide to materialographic specimen preparation.  
Contact your local dealer for a free copy of Metalog Guide™.

### **3. Trouble-Shooting**

Troubleshooting the AbraPlan-10 is achieved using popup messages displayed in the LCD display on the console unit. These messages are divided into several categories and are listed below with decreasing severity:

Fatal Errors

In the case of Fatal Errors, process cannot continue before an authorised technician has rectified the error. Turn off the unit at the main switch immediately. Do not attempt to operate the unit before a technician has rectified problem.

Errors

Errors must be rectified before process can continue.

Messages

Messages are intended to inform the operator of the machine's progress and advise about minor operational errors.

When a popup message is displayed, it will have the following format:

- a heading showing one of the categories listed above.
- an information field providing a message or instruction.

An example of a popup is shown below:

The information in the following table is divided by category, and then alphabetically by message.

*AbraPlan-10*  
*Instruction Manual*

<b>Error message</b>	<b>Explanation</b>	<b>Action required</b>
<b>Fatal error</b>		
15 V DC supply missing	No power to circuitry in AbraPlan-10	Contact Technical Support.
24 V DC supply missing	No power to circuitry in AbraPlan-10	Contact Technical Support.
No Serial communication	Communication problem between elements inside AbraPlan-10.	Contact Technical Support.
Program version mismatch	Incompatible software in AbraPlan-10.	Contact Technical Support.
<b>Error</b>		
Air pressure too low	There may be a leak in the hose or the compressor maybe defective.	Check the compressed air system for possible cause.
Disc Motor Overloaded	Load pressure is too high.	Reduce force value.
Dressing arm motor is stuck	This may be due to an internal electrical problem. Possible physical obstruction blocking arm from moving.	Fatal – contact Technical Support. Remove the obstruction.
Emergency Stop is active	The Emergency Stop is pushed in. This is displayed until action is taken.	Release the Emergency Stop.
Process has been stopped by lifting protective shield	The protective shield has been lifted.	Put the shield down into its operating position.
Pump Motor Overloaded	Pollutants in the coolant are clogging the pump.	Clean the cooling system and fill with fresh coolant mixture.
Removal rate too low, maximum permitted time has been exceeded	The machine has been grinding for 15 minutes without reaching the goal.	When using Stock Removal, check that the force was not altered during operation.
Sample holder unable to move upwards after process	There may be a problem with the compressed air system.	Check the compressed air system.
	There may be an internal electrical problem.	This is fatal – contact Technical Support.
Sample Motor Overloaded	Loading on the motor is too high.	Reduce the force value.

*AbraPlan-10*  
*Instruction Manual*

Message	Explanation	Action required
<b>Messages</b> (those in shaded boxes are safety precaution messages, where instructions are provided by the software)		
2-hand start released too early		Press the start buttons again, and hold them until the specimen holder begins to rotate.
Force regulation error	The specimen holder was moving when a loss in air pressure is experienced	Check the air pressure supply for leaks and faults.
	There may be an internal electrical problem.	This is fatal – contact Technical Support.
Grinding stone cover not closed	Message prior to dressing if the slide cover is not in position.	Move the slide cover all the way over to the left.
Grinding stone cover not open	Message received prior to grinding if slide cover is not in position.	Move the slide cover all the way over to the right.
Grinding stone must be replaced soon	Message received during the dressing process.	The grinding stone must be renewed as soon as possible.
Move dressing arm to top? Press F1 or F2	This message occurs during power up and after a new stone has been installed.	Move the dressing arm to the top to prevent damage to the new stone when dressing next time.
Process already stopping	Message if the Stop button is pushed when the process is already finished.	
Process in progress	Informs you that the process is active.	
Process is finished	Message at the end of the process.	
Process is stopping	This is recognition of the STOP command.	
Process has been emergency stopped, press F1 to raise sample mover	The process was stopped using the EMERGENCY STOP button. But the specimen holder is still down. Follow the instructions on screen to lift it.	Safety precaution message.  Follow the instructions on screen.
Process has been stopped by lifting the protective shield. Press F1 to raise sample mover, or lower shield again	The process was stopped by lifting the shield. But the specimen holder is still down. Follow the instructions on screen to lift it or continue.	Safety precaution message.  Follow the instructions on screen.
Shield not closed	The grinding process can not start because the shield is not down.	Lower the shield and start the process.
Start denied, select GRINDING Menu	The start buttons were pushed, but the process could not start from the current software menu.	Select the GRINDING Menu and try to start the process.
Stock removal sensor not installed	The software is set up for stock removal and the start buttons pushed, but the process won't start.	The stock removal option is not installed. Install it or select the timer option in the GRINDING Menu.
Water level in tank too low	More coolant is required – however the current process may be completed.	When the process is finished top up the tank with new coolant.

*AbraPlan-10*  
*Instruction Manual*

	<b>Explanation</b>	<b>Action</b>
<b>Physical Observations/Problems</b>		
No material removed	The grinding stone/diamond grinding disc is covered with material.	Dress the stone/disc. Replenish the cooling water volume.
	Insufficient grinding force.	Regulate the grinding force.
The specimen heats up	No cooling water.	Replenish the volume of cooling water in the recirculation unit.
		Cooling water pump is blocked.
Un-plane specimens	Stone not dressed.	Dress the stone. If this does not help, adjust the dressing arm (see the Maintenance section).
	Too few specimens in the specimen holder. Badly centred large specimen or specimen with too small a distribution in one direction.	Put blank specimen/s in the specimen holder.
Squeaking noise	V-belt slides.	Tighten the V-belt in accordance with the Maintenance section.
Hissing noise when the machine is in operation and/or is switched off	Leak in the air system.	Tighten the fittings and/or replace the defective air tubing.
Violent vibrations when the machine is running idle	The grinding stone is defective and out-of-balance.	Exchange the stone. Before doing so, try to turn the stone in relation to the turntable.
Machine very noisy when running idle	Axial bearing defective.	Exchange part.
	Spindle or motor bearings defective.	Dismantle the V-belt and start the motor. If the noise continues it means that the motor bearings should be changed.

## 4. Maintenance

### Daily Service

- Empty and clean the sieve in the Recirculation Cooling Unit.
- Refill the Recirculation Cooling Unit.

**IMPORTANT**

Always use the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to top up with Struers Additive each time you refill with water.

### Weekly Service

AbraPlan-10 and the Recirculation Cooling Unit should be cleaned regularly, in order to avoid damaging effects to the machine and the specimens from abrasive grains or metal particles.

When cleaning AbraPlan-10:

- Detach the cover over the grinding stone and remove any material that may have accumulated on the bottom of the drain pan.

### *Checking the Recirculation Cooling Unit*

The cooling unit should be checked for cooling water after 8 hours use or at least every week. The unit must be refilled if the flushing pump cannot reach the cooling water or if the cooling water is too contaminated. Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

To check the concentration of additive, use a refractometer. Concentration = Brix value. The concentration of additive should always lie between 2.7 and 3.3 %. Add Struers Additive for Cooling Fluid if the concentration is too low.

### Monthly Service

#### *Replacing the Cooling Water*

When cleaning the Recirculation Cooling Unit empty the tank of cooling water and clean the tank and the sedimentation container mounted under the discharge branch.

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.
- Place a disposable plastic insert (TREPO) in the tank and fold it over the edge.


**Please Note**

Too high a level of coolant in the tank might damage the pump.

- Fill the tank with 29.1 l of water and 900 ml Struers Additive.

## Adjusting the V-belts

The V-belts of the principal motor should be kept tight, in order to ensure transmission of the power to the grinding stone. Adjust the belts when necessary, or at least once every year. If a brief squeaking is heard when the grinding stone is started, the V-belts are too slack and should be tightened.

To adjust:

- Disconnect the main power.
- Remove the cover plate for the Recirculation Cooling Unit ①.
- Unscrew the 4 Allen screws securing the bottom cover and remove the cover.
- Unscrew the 2 Allen screws securing the cover plate and slide the plate off.
- Loosen the locking nut ②.
- Adjust the spindle ③ in such way that a load of about 15 N applied to the belt ④ at a point half-way between the two pulleys will cause the belt to move about 8 mm.
- Tighten locking nut after the adjustment.

## 5. Technical Data

Subject		Specifications	
		Metric/International	US
Grinding Stone/ Disc	Rotational speed	1450 rpm	1450 rpm
	Size	356 mm	14.0"
	Power consumption	4 kW	5.4 Hp
Specimens	Speed	150 rpm	150 rpm
	Direction	CCW	
	Force	50-700N	10-150 lbf
	Motor power consumption	0.37kW	0.5 Hp
Recirculation Cooling Unit	Power supply	3 (3L+PE)	3 (3L+PE)
	Power consumption	140W	0.2 Hp
	Width	500mm	19.7"
	Depth	400mm	15.8"
	Height	410mm	16.1"
	Weight	7.2kg	16.2lbs
Software and Electronics	LC Display	240x128 pixels	
	Controls	Touch pads/Push-turn knob	
	Memory	EPROM/RAM/NV-RAM	
Supply Voltage	Power consumption	4.4kW	
	No. of phases	3 (3L+PE)	
	Output, main motor	4.0kW	
	Voltage/frequency:		
	3 x 200 V / 50Hz	16.9A	
	3 x 200 - 210 V / 60Hz CSA	15.7A	
	3 x 220 - 230 V / 50 Hz	16.9A	
	3 x 220 - 240 V / 60 Hz	15.7A	
	3 x 380 - 415 V / 50 Hz	8.9A	
3 x 380 - 415 V / 60 Hz	10.3A		
3 x 460 - 480 V / 60 Hz CSA	8.5A		
Compressed Air	Compressed air supply	6-10bar	
Dimensions and Weight	Width	670mm	26.4"
	Depth	930mm	36.6"
	Height	1500mm	59.0"
	Weight	340kg	765lbs

*AbraPlan-10*  
*Instruction Manual*

<b>Subject</b>		<b>Specifications</b>	
		<b>Metric/International</b>	<b>US</b>
Environment	Safety standards	98/37/EØF EN292-1:1991; EN292-2:1991+A1; EN418:1992; EN574:1996; EN953:1997; EN954-1:1997; EN1050:1997; NFPA70:1996; UL508; NFPA79:1997 FCC part15	
	EMC	EN50081-1:1992; EN50082-1:1997; EN60204-1 1997	
	Noise level(idle)	77dbA	
	Surrounding temperature	5-40°C	
	Humidity	Max. 95%RH	

# Quick Reference Guide

## **Inserting the Specimen Holder**

- Position the specimen holder under the quick coupling.
- Press and hold the flange of the column down with the heel of your hand while guiding the pressure tap of the specimen holder into the coupling.
- Turn the specimen holder until the three pins engage with the corresponding holes.
- Release the flange.

## **Removing the Specimen Holder**

- Press upwards with your fingers to lift the specimen holder slightly. At the same time, press and hold down the flange with the heel of your hand.
- Using your fingers to support the specimen holder; lower it free of the coupling.
- Release the flange and completely remove the specimen holder.

## **Starting the Preparation Process**

- Insert the specimen holder.
- Ensure that the dressing arm is located in the “parked or resting” position on the left-hand side of the machine and that the sliding cover on the stone guard is located all the way to the right.
- Lower the shield.
- Enter the GRINDING menu and set/check the correct time and force.
- Start the grinding process.

## **Stopping the Preparation Process**

- When the time has elapsed, the grinding stone will automatically stop rotating and the specimen holder will return to its initial position.

## **Dressing the Grinding Stone**

- Press F1 to start the cooling water pump.
- Swing the diamond tool to a position over the grinding stone and move the diamond tool 2-6 times forwards and backwards across the grinding stone.
- Adjust the position of the diamond tool by pressing F3 or F4.
- When the dressing is complete, swing the diamond tool back to its “parked or resting” position.

## **Dressing the Diamond Grinding Disc**

- Mount 3 aluminium oxide dressing sticks in a sample holder and grind it for a few seconds.

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*AbraPlan-10*  
*Spare Parts and Diagrams*

---

**Always state *Serial No* and *Voltage/frequency*  
if you have technical questions or when ordering spare parts.**

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**Instruction Manuals:** Struers Instruction Manuals may only be used in connection with Struers equipment covered by the Instruction Manual.

**Service Manuals:** Struers Service Manuals may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

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#### **AbraPlan-10**

##### *Drawings*

AbraPlan-10, complete .....	15490001K
AbraPlan-10 Quick-release coupling, complete .....	15490007C
Quick-release coupling, complete.....	15090009B
Plate for electric connection, complete .....	15490016B
Stand, complete.....	15490010E
Drain assembly fitted to chassis .....	15490011E
Stone support disc.....	15490005D
Stone cover .....	15490008A
Flushing hose, complete .....	15490009B
Plate with process print and pneumatic, complete.....	15490031E
Front plate, complete.....	15490060C
Air filter/regulator unit, complete.....	15490015B
Contacting box, complete .....	15490033D
Cabinet for drain tray, complete .....	15490050F
Column guide, complete .....	15490020G
M1 Main motor, complete.....	15490021C
Bearing unit sub assembly.....	15090040A
M2 Sample Mover motor, complete .....	15490045B
Dresser & lift cylinder, complete .....	15490070F
Dresser arm.....	15490072F
Stock removal unit.....	15490083D
TREJA, Rolling pallet for cooling unit, low version.....	15490080

##### *Diagrams*

Block diagram .....	15493050A
Circuit diagram.....	15493100D
Control voltage diagram (3 pages) .....	15493120F
Functional diagram .....	15493130C
Air diagram .....	15492000B

Some of the drawings may contain position numbers  
not used in connection with this manual.

*AbraPlan-10*  
*Spare Parts and Diagrams*

**Spare Part list for AbraPlan-10**

**Drawing**  
 15490001

<b>Pos.</b>		<b>Cat no.</b>
	<b>AbraPlan-10, complete</b>	
20	Drain Tray fitted to Cabinet	15490011
70	Column guide, complete	15490020
110	Rubber bushing	15090690
170	Sample mover motor cover plate	15490306
180	M2 Sample Mover motor, complete	15490045
183	Safety switch 1NO+1NC (4R)	RSS00008
190	2.00A Fuse Midget 10x38	2FC10020
192	4.00A Fuse Midget 10x38	2FC10040
195	1.00A Fuse Midget 10x38	2FC10010
220	Process Control PCB + Pneumatic subassembly.	15490031
250	Plate for electrical connection, complete	15490016
280	Front panel, complete	15490060
425	Lifting iron bar with screws AbraPlan-10	R5494100
470	Quick-release coupling, complete	15490007
474-477	Screw & key for quick-release coupling	15092920
630	Flushing gun, complete	15490009
640	Spiral hose, ø32	2NU29401
650	Hose Clamp ø25-40	2NS22540
660	Dresser assembly, complete	15490070
740	Safety switch 1NO+1NC (4R)	RSS00008
900	Hose nipple	13688045
910	Allen key 10	2GR00100
920	Fork wrench NV13	2GR00213
930	Fork wrench NV24 DIN 894	2GR00224
940	Safety glasses	2GU10318
950	TRAJE Rolling Pallet for Cooling Unit, low version	05496901
960	ABPEU Stock removal unit	05496902
970	PAMWA water level sensor	03756912
980	TRECI Cooling unit 3X415-480V/50-60Hz	02536154
990	TRECI Cooling unit 3X345-415V/50-60Hz	02536149
1000	TRECI Cooling unit 3X200-240V/50-60Hz	02536135

15490007

	<b>AbraPlan-10</b>	
	<b>Quick-release coupling, complete</b>	
10	Quick release coupling, complete	R5090007
30-50	Splash ring w. screws	260MP060

*AbraPlan-10*  
*Spare Parts and Diagrams*

**Spare Part list for AbraPlan-10**

<b>Drawing</b>	<b>Pos.</b>	<b>Cat no.</b>
15090009	<b>Quick-release coupling, complete</b>	
	5	Quick release coupling, complete R5090007
	6	Screw & key for quick-release coupling, complete 15092920
	30	Pressure spring Ø4.3X0.7 12600718
	40	Driving pin 15490420
	60	Key /fitted key, type A 12600723
	65	Ball KU 5.556 (ø7/32 in) RS 2BA00055
15490016	<b>Plate for electric connection, complete</b>	
	10	Plate for electric connection, complete R5090015
	215	Terminal strip SAK 6/35 PA 2XK41012
	216	Ground terminal 2XR42617
15490010	<b>Stand, complete</b>	
	20	Machine leg M75-16-070 2GB40075
	150	Air filter/regulator unit, complete 15090032
15490011	<b>Drain assembly fitted to chassis</b>	
	80	Drain tray 15490057
	90	Stone support disc 15490005
	105	Upper flange for grinding stone R4590026
	120	Stone cover assembly 15490008
	130	Knurled handle, M8 DIN466-NI 2GH00209
	160	Safety switch 1xNO+1xNC RSS00008
	220	Elbow 87 for hose ø38/ pipe socket ø40 2NG20587
	230	Hose clamp 40-60 mm 2NS24060
240	Drain hose ø51 2NU30251	
15490005	<b>Stone support disc</b>	
	10	Stone support disc 15490005
15490008	<b>Stone cover</b>	
	80	O-RING 12.42-1.78 72 NBR 872 2IO17817
	110	INA-Sealing ring G 10X17X3 2II01017

*AbraPlan-10*  
*Spare Parts and Diagrams*

**Spare Part list for AbraPlan-10**

<b>Drawing</b>	<b>Pos.</b>	<b>Cat no.</b>	
15490031	<b>Plate with process print and pneumatic, complete</b>		
	20	PCB A2 AbraPlan+I458-10 15493002	
	21	Set of Proms, latest version 15492900	
	22	2.5AT Fuse Sub-Miniature 2FU32900	
	23	4.0AT Fuse Sub-Miniature 2FU32999	
	60	Pressure governor, 5-8.5 bar 1/4in 2YR00001	
	70	Gasket, PVC O-1/8 2IF00011	
	100	Gasket, PVC O-1/4 2IF00012	
	150	3/2 solenoid valve 24V DC NORM.Å. 2YM10120	
	190	Gasket, PVC M5 2IF00010	
	220	3/2 solenoid valve 24VDC 2YM10124	
	240	Angle for quick-coupling, rotatable ø5-M5 2NF10081	
	270	Gasket, PVC M5 2IF00010	
	280	Clamping ring 14590110	
	285	PU hose PU-3 blue ø2.9/ø4.3 2NU14031	
	290	Nylon hose, SUPERFLEX ø5/ø3.2 2NU12445	
	310	Neoprene nipple ø36/ø47/ø54-2.5 2GK90457	
	330	20 Ohm Break resistor HSD40A 2RK05602	
	15490060	<b>Front plate, complete</b>	
		10	Front plate with foil and window pane 15492909
60		Display with plug 14930045	
80		PCB A1 control panel AbraPlan-10 15093001	
81		Set of Proms, latest version 15492900	
110		Light diode PCB 15022396	
150		2 channel opt. encoder w. pressure 24p 2HR12411	
160	Turn/push-button R5090600		
1540032	<b>Air filter/regulator unit, complete</b>		
	40	Air filter, air regulation EAW3000-F02D-6 2YF00005	
	90	Gasket, PVC 1/4" 2IF00012	
	110	Gasket, PVC 1/8" 2IF00011	
	140	PVC-Hose, clear 13/32"-Ø10 2NU19313	
150	Nylon hose, SUPERFLEX ø5/ø3.2 2NU12445		

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

CONNECTION TO SAFETY TRANSFORMER 2MT72034

IEC742 SAFETY TRANSFORMER

SUPPLY OUT: 2 X 24Vac, 4AT

Note: Refer to drawing 15493100-dsn for correct wiring.

M1 3 PH Sample mover motor 370W / 300rpm

M2 3 PH Disc/grind motor 4000W / 2850rpm

M3 3 PH Recirculation pump 160W

PCB 12050400

MACHINE CONTROL-BOARD, A2.

PCB 120504XX

CPU-BOARD, A1.


SUPPLY IN: 2 X 24Vac, 4AT

- Inputs:**
- Status twohand relay
  - Emergency stop
  - Temperature disc motor
  - Temperature sample mover motor
  - Temperature recirk. motor
  - Level sensor recirk.
  - Sample mover bottom position
  - Sample mover top position
  - Dresser in parked position
  - Lid open
  - Lid closed
  - Visor open
  - Removal measurementunit installed
  - Absolute pressure button
  - Absolute pressure top
  - Voltage supply 1 (input)
  - Voltage supply 2 (+15V)
  - Voltage supply 3 (+24V)
  - Removal measurement
  - Motor current
  - Dresser removal
  - Speed sample mover motor (freq. conv.)
  - Speed disc motor (freq. conv.)

- Outputs:**
- 7 x contactors:
    - Main
    - Sample mover 150 rpm
    - Disc motor 1450 rpm
    - Disc motor star
    - Disc motor delta
    - DC brake disc motor
  - Recirculation
  - 5 x 24VDC air valve:
    - Sample mover up
    - Air outlet 1
    - Air outlet 2
    - Sample mover down
    - Dresser up
  - 1 x 24VDC dresser motor



OPTION: Freq. converter.

STRUERS A/S VALHOESALLE 178 DK-2610 ROEDVRE DENMARK PHONE: + 45 3670 3500			
ABRAPLAN 10, COMPLETE BLOCK DIAGRAM			
FILE NAME: 549-3050.DSN PAGE1.DSN	Size A2	CAGE Code <Cage Code>	DWG NO 549-3050
Scale	Scale	Sheet	Rev A
Wednesday, January 10, 2001	Scale	SLN / SLN	Sheet 1 of 1


MUST BE EXTERNALLY FUSED  
EXTERNAL FUSES = 25AT

MOTOR COUPLING FOR M2

ALL OTHERS  
AS DELIVERED



ATTENTION!  
Do not use any shorting links

COLOR CODE:  
BK = black  
BN = brown  
RD = red  
OG = orange  
YE = yellow  
GN = green  
BU = blue  
VT = violet  
GY = grey  
WH = white  
PK = pink  
BE = beige  
RO = rose

INPUT VOTAGE	F1+F2+F3	F4+F5	F6+F7+F8	W1	W2	W25
3 x 200V / 50Hz	3 x 2AT	2 x 4AT	3 x 1AT	2,5mm <sup>2</sup>	1,5mm <sup>2</sup>	1,5mm <sup>2</sup>
3 x 200-210V / 60Hz	3 x 2AT	2 x 4AT	3 x 1AT	AWG = 12	AWG = 14	AWG = 14
3 x 220-230V / 50Hz	3 x 2AT	2 x 4AT	3 x 1AT	2,5mm <sup>2</sup>	1,5mm <sup>2</sup>	1,5mm <sup>2</sup>
3 x 220-240V / 60Hz	3 x 2AT	2 x 4AT	3 x 1AT	AWG = 12	AWG = 14	AWG = 14
3 x 380-415V / 50Hz	3 x 1AT	2 x 2AT	3 x 1AT	2,5mm <sup>2</sup>	1,5mm <sup>2</sup>	1,5mm <sup>2</sup>
3 x 380-415V / 60Hz	3 x 1AT	2 x 2AT	3 x 1AT	AWG = 12	AWG = 14	AWG = 14
3 x 460-480V / 60Hz	3 x 1AT	2 x 2AT	3 x 1AT	AWG = 12	AWG = 14	AWG = 14


Rev. A: Baan Rev. B: Wiring size changed Rev. C: Coupling on M1 changed	STRUERS A/S VALHOEJS ALLÉ 176 DK-2610 ROEDOVRE DENMARK PHONE: +45 3870 3500
<b>ABRAPLAN 10, CIRCUIT DIAGRAM MAIN VOLTAGE, MAIN SUPPLY CIRCUIT.</b>	
FILE NAME: 549-3100_DSN PAGE: 1 OF 1	Size A2 CAGE Code DWG NO <b>549-3100</b> IRev
Monday, June 17, 2002	Scale SLN / SLN Sheet 1 of 1


COLOR CODE:  
 BK = black  
 BN = brown  
 RD = red  
 OG = orange  
 YE = yellow  
 GN = green  
 BU = blue  
 VT = violet  
 GY = grey  
 WH = white  
 PK = pink  
 BE = beige  
 RO = rose

Rev. A: New Rev. B: Wiring changed Rev. C: Wire changed for dresser PCD Rev. D: Wire changed for removal design 1344 Rev. E: Wire Color to V2 & V3 (wiring changed) (2006/12/05) Rev. F: Wire changed to V2 (2006/12/05)	STREIBER AS VAL-VE-DE-ALLE 116 DK-2010-ROD-CO-VE 2006/04/04 PDK-NL-45-3070-3500	ABRAPLAN 10, CONTROL VOLTAGE DIAGRAM MACHINE CONTROL-BOARD		Rev. No.
F111: 04/01: 24-1-2010: 008	Rev. A1	Rev. Code	DWG. NO.	549-3120
Friday, May 29, 2010	Scale	0:1.00 / 0:1.00	Sheet	2 of 3

COLOR CODE:  
 BK = black  
 BN = brown  
 RD = red  
 OG = orange  
 YE = yellow  
 GN = green  
 BU = blue  
 VT = violet  
 GY = grey  
 WH = white  
 PK = pink  
 BE = beige  
 RO = rose


Rev. A: New Rev. B: Wiring changed Rev. C: Wire changed for dresser PCB Rev. D: Wire changed for removal design (S10) Rev. E: Wire Color to V2 & V3 (wiring changed) (S10) Rev. F: Wire changed to V2 (S10)	STREIBER AS VALVE-REIN. ALLE 116 DK-2010 (ROD) (V2) 20050404 PDK-NL_45_2010_3500	<b>ABRAPLAN 10, CONTROL VOLTAGE DIAGRAM MACHINE CONTROL-BOARD</b>		DWG NO <b>549-3120</b>	Proj <b>F</b>
FILE NAME: 1_24-1-2010_008 DATE: 01-24-2010	SHEET <b>A1</b>	SCALE <b>1:1</b>	DWG NO <b>549-3120</b>	Date <b>Friday, May 29, 2010</b>	Page <b>2 of 2</b>

COLOR CODE:  
 BK = black  
 BN = brown  
 RD = red  
 OG = orange  
 YE = yellow  
 GN = green  
 BU = blue  
 VT = violet  
 GY = grey  
 WH = white  
 PK = pink  
 BE = beige  
 RO = rose

Rev. A: Baan Rev. B: Wiring changed Rev. C: Connection between L1 and L2 removed at 380-420V and 460-480V	STRUERS A/S VALHOEJS ALLÉ 176 DK-2610 ROEDOVRE DENMARK PHONE: + 45 3670 3500
ABRAPLAN 10, FUNCTIONAL DIAGRAM, SAFETY CIRCUIT	
FILE NAME: 549-3130_DSN PAGE: 8/8	Size A2 CAGE Code DWG NO 549-3130 Rev I
Thursday, April 24, 2003	Scale SUN / KFC Sheet 1 of 1

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Herewith declares that

<i>Product Name:</i>	AbraPlan-10
<i>Type No:</i>	549
<i>Machine Type:</i>	Grinding machine

is in conformity with the provisions of the following directives:

**Safety of Machinery** 98/37/ EEC according to the following standard(s):  
EN ISO 12100-1:2003, EN ISO 12100-2:2003, EN ISO 13850:2006, EN574:1996; EN953:1997;  
EN954-1:1997; EN1050:1997; EN60204-1:1997.

**EMC-Directive** 2004/108/EEC according to the following standard(s):  
EN61000-6-1:2002, EN61000-6-3:2002.

**Low Voltage Directive** 2006/95/EEC according to the following standard(s):  
EN60204-1:1997.

**Supplementary Information** The equipment complies with the American standards:  
FCC part 15, subpart B og UL508 ; NFPA 70:1996; NFPA 79:1997

**The above has been declared according to the global method, module A**

Date: 29.11.2007

Christian Skjold Heyde,  
Vice President, R& D and Production, Struers A/S

**Hersteller** Struers A/S  
Pederstrupvej 84  
DK-2750 Ballerup, Danmark  
Telefon +45 44 600 800

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<i>Produktname:</i>	AbraPlan-10
<i>Typennr.:</i>	549
<i>Maschinenart:</i>	Schleifmaschine

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# AbraPlan-10

 **Struers**

**An automatic,  
high-capacity machine for  
fast and efficient plane  
grinding**

**AbraPlan-10 is a powerful and fast plane grinding machine ideal for laboratories or workplaces carrying out quality control testing of very large specimens or a high volume of specimens**

**Unique advantages:**

- Fast and efficient plane grinding providing high planeness and reproducibility, thus saving valuable time on the subsequent preparation steps
- Easy to operate without requiring extensive training
- Sturdy and powerful design ensuring troublefree operation year after year
- Semi-automatic dressing arm monitoring the wear on the grinding stone
- High level of operator safety
- Built-in flush hose for easy maintenance
- Controlled material removal (Optional)

**The AbraPlan-10** is an ideal plane grinding machine for users with very large specimens or with a high volume of specimens.

On the 356 mm / 14" diameter grinding stones, 200 mm / 8" diameter sample holders can be utilized efficiently taking up to 12 specimens at a time.

Whenever capacity or preparation time is crucial, AbraPlan-10 is the obvious choice. One AbraPlan-10 can easily provide enough work for two or three polishing machines or be employed as a pre-grinder for the fully-automatic MAPS system.


### Fast and Efficient Plane Grinding

The grinding stone rotates at 1450 rpm guaranteeing a very high material removal. Using  $Al_2O_3$  and SiC grinding stones, AbraPlan-10 can grind both ferrous and non-ferrous metals. For ceramics and sintered carbides a resin bonded diamond grinding disc is available.

With typical grinding times of 20-30 seconds, a high quality surface finish is achieved quickly. Specimen planeness and reproducibility is unmatched on the materialographic market, saving time on the subsequent fine grinding and polishing steps.

### Easy to Operate

The AbraPlan-10 is easy to operate and extensive training is therefore unnecessary. A large display offers the operator a good overview of the process parameters and the push/turn knob makes it easy to select and edit parameters.

*AbraPlan-10 is the strongest plane grinding machine available on the materialographic market*

*The operator can monitor the thickness of the stone through its entire lifetime*


A countdown feature allows the operator to see the remaining grinding time or remaining material removal.

### Sturdy and Powerful Design

AbraPol-10 is the strongest plane grinding machine available on the materialographic market. The powerful motors, the sturdy design and the wellknown Struers technology contribute to many years of trouble-free use.

### Optional Material Removal Unit

In some situations grinding for a specified grinding time is not a good


*A large display offers a good overview of the process parameters*

idea. With small samples too much material might be removed, while large samples or hard materials may barely have been touched. Therefore, a sensor to accurately measure the material removal has been developed. The Accurate Material Removal Unit allows you to preset the amount of material that has to be removed. The removal can be set from 0.01 mm to 2.00 mm in increments of 10  $\mu$ m.

### Semi-Automatic Dressing

The AbraPlan-10 is provided with a semi-automatic diamond dresser that keeps the grinding stone plane and sharp. The position of the dressing arm and the incremental dressing distance are software controlled. This means that the operator can monitor the thickness of the stone through

*Inside AbraPlan-10 there is room for a recirculation cooling unit*

its entire lifetime. The up and down movement of the dressing arm is pneumatically supported. The operator will therefore not have to apply any force, but can simply swing the dressing arm with the diamond tool across the grinding stone.

**High Safety Level**

Straus recognizes the importance of a safe working environment and has equipped AbraPlan-10 with several safety features.


The grinding stone and the sample holder are totally encapsulated, ensuring the operator cannot get his hands caught in the holder or by the stone during grinding or dressing. A shield offers additional operator protection. The AbraPlan-10 cannot be started unless these protectors are closed correctly, and the machine will stop immediately if an attempt is made to open them during operation. AbraPlan-10 is provided with a two-hand start so the operator's hands

are outside the machine when starting, and an emergency stop is located at the front of the machine.

**Built-in Flush Hose**

The flushing unit allows for easy cleaning of the grinding area, samples and the grinding chamber. A recirculation cooling unit can be connected for cost-efficient water utilization and to collect removed abrasives.

*One AbraPlan-10 can easily provide enough work for two or three grinding and polishing machines such as the automatic AbraPol-10, a RotoSystem or MAPS*

## Technical Data

<b>Grinding wheel</b>	Rotational speed Size	1450 rpm 356 mm / 14"
<b>Specimen holder</b>	Speed Force Max. size sample holder	150 rpm 50-700 N 200 mm / 8"
<b>Software and Electronics</b>	LC Display Controls	240 x 128 characters Touch pad and Push-turn knob
<b>Environment</b>	Noise level Surrounding temperature Humidity	70 dB(A) 5-40 C Max. 95 % RH
<b>Supply voltage</b>	Power consumption No. of phases Output, main motor	4.5 kW 3 (3L+PE) 4 kW
<b>Compressed air</b>	Compressed air supply	6-10 bar
<b>Dimensions and weight</b>	Width Depth Height Weight	670 mm / 26.5" 939 mm / 37" 1500 mm / 59" 340 kg / 750 lbs

*Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)*

## Specifications

### AbraPlan-10

Automatic plane grinding machine for specimens in specimen holders 160 mm / 6.3" and 200 mm / 8" diameter. With diamond dresser. Recirculation cooling unit (TREC1), grinding stone and specimen holders are ordered separately.

**Code**

ABPLA

## Accessories

### Accurate Material Removal Unit

Stops grinding process after a specified amount of material has been removed.

ABPEU

### Recirculation Cooling Unit

For water cooling on various equipment.

Capacity 30 l.

TREC1

### Rolling Pallet for TREC1

TREJA

### Water Level Indicator

Indicating the level of water in the recirculation cooling unit TREC1.

PAMWA

### Grinding Stone

Grit 60. Aluminium oxide grinding stone for plane grinding of soft and ductile metals (<HV 250). 356 mm / 14" diameter.

ABGOF

### Grinding Stone

Grit 150. Aluminium oxide grinding stone for plane grinding of harder metals (>HV 250). 356 mm / 14" diameter.

ABGAL

### Grinding Stone

Grit 150. SiC grinding stone for plane grinding of non-ferrous metals. 356 mm / 14" diameter

ABGAS

### Diamond Grinding Disc

Grit 120. Resin bonded disc, grooved in a square pattern, for grinding of ceramics and sintered carbides. 356 mm / 14" diameter.

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*Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.*

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