# ProntoPress-10

## Instruction Manual

Manual No.: 15067003

Date of Release 11.01.2008

*ProntoPress-10*  
*Instruction Manual*

<b>Table of Contents</b>	<b>Page</b>
User's Guide .....	1
Reference Guide.....	21
Quick Reference Card .....	54

---

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

**Instruction Manuals:** Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

**Service Manuals:** Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2008.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup  
Denmark  
Telephone +45 44 600 800  
Fax +45 44 600 801

---

## **ProntoPress-10**

### **Safety Precaution Sheet**

#### **To be read carefully before use**

1. The operator should be fully aware of the use of the machine according to the Instruction Manual.
2. The machine must be placed in a well ventilated room on a working table with adequate height for handy operation.
3. Be sure that the actual voltage corresponds to the voltage stated on the back of the machine and on the heating/cooling unit. The machine must be earthed.
4. Be sure that the water connections are mounted correctly and without leaks. The main water supply should be turned on when the machine is in use.  
Struers recommend that the mains water supply is shut off or disconnected if the machine is to be left unattended.
5. Be sure that the outlet hose is attached to the water outlet system in a safe way.
6. Be sure that the mounting unit is correctly assembled on the press before starting the process.
7. Be sure that the top closure with upper ram is either correctly mounted on the mounting cylinder or completely removed from the mounting cylinder before starting the press.
8. Do not operate the mounting press with a higher force/pressure than recommended for the actual cylinder diameter and resin in Struers Application Guide for Hot Mounting.
9. Following a heating cycle, ensure the mounting cylinder is cooled for a minimum of two minutes before opening.
10. Do not operate the machine whilst assembling or disassembling the mounting unit.

---

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

---

# User's Guide

Table of Contents	Page
<b>1. Getting Started</b>	
Checking the Contents of Packing.....	3
ProntoPress-10.....	3
Mounting Unit.....	3
Unpacking ProntoPress-10.....	3
Placing ProntoPress-10.....	3
Transport Screw Cap.....	3
Getting Acquainted with ProntoPress-10.....	4
Noise Level.....	5
Supplying Power.....	5
Supplying Water.....	5
Water Inlet.....	5
Water Outlet.....	5
Assembling the Mounting Unit.....	6
Removing the Cover.....	6
Installing the lower ram.....	6
Installing the Mounting Unit.....	6
Installing the Cover.....	7
Installing the Dust Protection Ring.....	7
Installing the Swivel Arm.....	7
Mounting of the Top Closure.....	7
Disassembling the Mounting Unit.....	7
Removing the Swivel Arm.....	7
Descending the Lower Ram.....	7
Removing the Dust Protection Ring.....	7
Removing the Cover.....	8
Removing the Mounting Unit.....	8
Removing the lower ram.....	8
Changing the Mounting Unit.....	8
Software Settings.....	9
Configuration Menu.....	9
Setting the Language.....	10
Setting the Cylinder Diameter.....	11

## **2. Basic Operations**

Using the Controls.....	12
Front Panel Controls of ProntoPress-10 .....	12
Groups of Keys.....	12
Acoustic Signals.....	12
Display .....	14
Main Switch .....	14
Mounting Method Options .....	15
Select a Mounting Method .....	15
Method Edit .....	17
Placing the Specimen .....	18
Pouring Resin over the Specimen.....	18
Placing Two Specimens .....	19
Installing the Top Closure.....	19
Starting the Mounting Process .....	20
Display during the Mounting Process .....	20
Stopping the Mounting Process .....	20
Removing the Top Closure .....	20

## 1. Getting Started

### Checking the Contents of Packing

#### *ProntoPress-10*

In the packing box you should find the following parts:

- 1 ProntoPress-10 machine
- 1 Swivel arm
- 1 Pressure hose
- 1 Filter gasket
- 1 Gasket
- 1 Reduction ring with gasket
- 2 Measuring spoons for mounting resin
- 1 Funnel
- 1 Air filter
- 1 Set of Instruction manuals

#### *Mounting Unit*

- 1 Mounting unit
- 1 Top closure with upper ram
- 1 Lower ram
- 1 Piston pin
- 1 Dust protection ring
- 1 Mould release agent (FASTI)
- 1 Scraper (PROAN)


### Unpacking ProntoPress-10

ProntoPress-10 is detached from the bottom of the packing case by removing the four screws from below.

### Placing ProntoPress-10

ProntoPress-10 should be placed on a steady table with an adequate working height. The machine must be placed close to the power supply, water mains and water outlet facilities. If water recirculation is used, there must be room under the table for the Recirculation Cooling Unit (TRECI).

### Transport Screw Cap


Carry out the following procedure before using the machine for the first time Exchange the transport screw cap ① with the enclosed air filter ②, to equalise the pressure in the hydraulic system.

The transport screw cap ① is situated underneath the cover for the mounting unit.

**Getting Acquainted with  
ProntoPress-10**

Take a moment to familiarise yourself with the location and names of the ProntoPress-10 components.

- ① Front panel
- ② Top closure
- ③ Main switch
- ④ Mounting unit
- ⑤ Dust protection ring
- ⑥ Cover for mounting unit
- ⑦ Swivel arm for top closure

## Noise Level

The noise level of the machine is 60 dB (A) measured when the pump is running, at a distance of 1.0 m/39.4" from the machine.

## Supplying Power

ProntoPress-10 is factory mounted with an electric cable. Mount a plug on the cable:

- Brown: phase
- Blue: neutral
- Yellow/green: earth


### **IMPORTANT**

Check that the mains voltage corresponds to the voltage stated on the type plate on the machine. Do not operate the machine before the transit screw cap has been exchanged with the enclosed air filter.

## Supplying Water

Water for cooling may be supplied from the water mains or a Recirculation Cooling Unit.

### *Water Inlet*

Mount the pressure hose onto the water inlet tube (pos.1) on the back of the machine:

- Insert the filter gasket in the coupling nut with the flat side against the pressure hose.
- Tighten the coupling nut completely.

- ① Water inlet
- ② Water outlet tube

Mount the other end of the pressure hose on the water mains tap for cold water:

- Mount the reduction ring with gasket on the water mains tap, if necessary.
- Introduce the gasket and tighten the coupling nut completely.

### **IMPORTANT**

Only connect to cold water


### *Water Outlet*

- Lead the water outlet tube (pos. 2) to the drain and be absolutely sure to place the hose with a steady downward slope and without any obstructions.
- Attach the water outlet tube to the water outlet system.

## Assembling the Mounting Unit

### Removing the Cover

- Loosen the catch (pos. 4).
- Press gently on both sides of the cover, to release it, and lift it up.

④ Catch


### Installing the lower ram

- Place the lower ram on the top of the piston rod.
- Turn the lower ram so that the holes in the axle journal on top of the rod and ram line up, and insert the piston pin.
- Make sure that the ends of the pin do not protrude.

#### **Important**

The axle journal is part of the machine's safety mechanism to protect against damage to the machine. If it breaks, it must only be replaced with a Struers replacement part or the safety mechanism may not function.


### Installing the Mounting Unit


- ① Mounting cylinder
- ② Fixation screw
- ③ Upper quick coupling for cooling water
- ④ Lower quick coupling for cooling water
- ⑤ Plug

#### **IMPORTANT**

Do not operate the machine whilst installing the mounting unit

- Unscrew the fixation screw about 10 mm.
- Place the mounting cylinder over the lower ram with the fixation screw in the position shown.
- Turn the mounting unit in a clockwise direction until it stops.
- Tighten the fixation screw completely.
- Mount the tube with the straight quick coupling on the lower quick coupling of the mounting unit (pos. 4) and push to connect. Make sure that the ring comes to stop in the extreme end of the quick coupling.
- Mount the tube with the elbow quick coupling on the upper quick coupling of the mounting unit (pos. 3) and push to connect. Make sure that the ring comes to stop in the extreme end of the quick coupling.
- Mount the plug (pos. 5) in the socket. Tighten the coupling nut.

④ Catch

*Installing the Cover*

- Introduce the barb in the front of the opening for the cover.
- Introduce the barbs on both sides of the cover in the cabinet. Press the sides of the cabinet gently.
- Press the rear end of the cover in so that the catch engages.

*Installing the Dust Protection Ring*

Place the dust protection ring around the mounting cylinder. The concave side should face upwards.

*Installing the Swivel Arm*

Insert the rod in the hole at the top of the cover for the mounting unit.

*Mounting of the Top Closure*

Introduce the top closure in the hole at the top of the swivel arm.

**Disassembling the Mounting Unit**

*Removing the Swivel Arm*

Lift the swivel arm out of the hole at the top of the cover for the mounting unit.

*Descending the Lower Ram*


- Switch on the main power of the machine.
- Press the key RAM DOWN ▼ to lower the ram to its lower limit.

*Removing the Dust Protection Ring*

Lift up the dust protection ring.


### Removing the Cover

- Press the catch (pos. 4).
- Press gently on both sides of the cover, to release it, and lift it up.

④ Catch

### Removing the Mounting Unit

- ① Mounting cylinder
- ② Fixation screw
- ③ Upper quick coupling for cooling water
- ④ Lower quick coupling for cooling water
- ⑤ Plug

**IMPORTANT**

Do not operate the machine whilst removing the mounting unit

- Disconnect the plug (pos. 5) in the socket. Loosen the coupling nut and pull.
- Disconnect the tube with elbow quick coupling on the upper quick coupling on the mounting unit (pos. 3). Pull the ring. Wait 5 seconds to allow the water to flow out of the cooling coil.
- Disconnect the tube with straight quick coupling on the lower quick coupling on the mounting unit (pos. 4). Pull the ring.
- Unscrew the fixation screw about 10 mm.
- Turn the mounting unit in an anti-clockwise direction until it stops.
- Lift the mounting unit.

### Removing the lower ram

- Push the piston pin out of the lower ram.

**Important**

The axle journal is part of the machine's safety mechanism to protect against damage to the machine. Please see: Maintenance. If it breaks, it must only be replaced with a Struers replacement part or the safety mechanism may not function.

### Changing the Mounting Unit

Follow the instructions for Removing the Mounting Unit and Installing the Mounting Unit.

**Software Settings**  
*Configuration Menu*

Switch on the power at the main switch. The following display will appear briefly:

```
ProntoPress-10  
Version: 2.0
```

The mounting press will start in the same state, as when the voltage was last switched off. The display to appear could be the Select Method Group Menu. If the heading in the display is different, press Esc until the Select Method Group Menu appears.

The Select Method Group Menu is the highest level in the menu structure. From here you can go to configuration, manual functions, select mounting methods in the database and open the mounting methods.

```
SELECT METHOD GROUP ↓  
1. STRUERS                2. GROUP NUMBER 2  
3. GROUP NUMBER 3          4. GROUP NUMBER 4  
F1:CONFIG.      F3:MAN. FUNC.
```

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

Setting the Language

Press MENU UP/DOWN ▲ ▼ to select the Language parameter.

Press ENTER to activate the Language Options.

Press MENU UP/DOWN ▲ ▼ to select the Language you prefer.

Press ENTER to accept a language. The Configuration Menu now appears in the language you have chosen.

**PLEASE NOTE**  
Voltage versions for 220 - 240 V:  
include the languages English, French and German  
Voltage versions for 100 - 120 V:  
include the languages English, French and Japanese

### Setting the Cylinder Diameter

The cylinder diameter has to be set, in order to have the right values of force or pressure. When a cylinder diameter is changed, the database values of force or pressure in the methods will automatically change to maintain the conditions in the mounting cylinder.

**Note:** The counter, which warns the user to clean the rams after 200 mountings, will be reset when a new cylinder diameter is selected. It is therefore important to clean away any old cured resin beneath the lower ram when changing cylinder diameter. Please see the instructions in the 'Maintenance' section.

 Press MENU UP/DOWN  to select the Cylinder parameters.

↓

CONFIGURATION		↓
Cylinder Diameter:	50 mm	
Temperature Unit:	CELSIUS	
Force or Pressure:	FORCE	Unit: kN


↓

 Press ENTER  to activate the Cylinder Diameter Menu.

↓

CYLINDER DIAMETER			
25 mm	30 mm	40 mm	50 mm
1 1/4"	1 1/2"		

↓

 Press MENU UP/DOWN  to select the diameter for the cylinder, which is mounted.

↓

 Press ENTER  to accept the diameter.

↓


 Press Esc to return to the Select Method Group Menu.

Setting the Temperature Unit, Force or Pressure Mode, Force Unit, Pressure Unit and Countdown Mode: See the section Advanced Operations.

## 2. Basic Operations

### Using the Controls

Front Panel Controls of  
ProntoPress-10

### Groups of Keys









- A** Programming and monitoring functions.
- B** Operation controls for the mounting cylinder.

### Acoustic Signals

-  *Short Beep:* when a key is pressed, a short beep indicates that the command has been accepted.
-  *Long Beep:* a long beep indicates that the key is inactive at that moment.
-  *A sequence of four double beeps* indicates an error.
-  *Three long beeps:* indicate that the mounting process is finished.


*ProntoPress-10*  
*Instruction Manual*

Pos. No.	Key	Function	Pos. No.	Key	Function
① MAIN SWITCH		Turns the main power of the machine on/off. The main switch is located at the right hand side of the cabinet.	⑫ FUNCTION KEY		Control for various purposes.
② POWER		Lights up when the main power is on.	⑬ FUNCTION KEY		Control for various purposes.
③ START		Starts the automatic mounting process.	⑭ FUNCTION KEY		Control for various purposes.
④ STOP		Stops the mounting process. Stops upward or downward movement of the lower ram.	⑮ MENU UP		Moves the cursor/page up or increases the parameter value.
⑤ RAM UP		Starts the upward movement of the lower ram. The ram automatically stops when the upper limit is reached.	⑯ MENU DOWN		Moves the cursor/page down or decreases the parameter value.
⑥ RAM DOWN		Starts the downward movement of the lower ram. The ram automatically stops when the lower limit is reached.	⑰ Esc		Leaves the present menu, moves back one step or cancels changes.
⑪ FUNCTION KEY		Control for various purposes.	⑱ ENTER		Selects a menu or allows editing of the marked parameters.

*ProntoPress-10*  
*Instruction Manual*


*Display*

- A** Heading.
- B** Flashing text. Cursor.
- C** Function key options.
- D** Arrow indicates that there are more lines in the picture.

**Please Note**  
The examples of display screens in this Instruction Manual show a number of possible texts. The actual display screen may differ from the examples in the Instruction Manual.

*Main Switch*

## Mounting Method Options

The mounting press will start in the same state, as when the voltage was last switched off. If a mounting method was selected at that time, this method will appear on the display.

```
Grp: 2 Method: 3 Resin:MultiFast
HEAT: Time: 7:00 T: 180°C F: 20 kN
COOLING: Rate: HIGH Time: 3:00
F1:SELECT PREHEATING
```

- If this mounting method can be used again, go to the section Starting the Mounting Process.
- If a new mounting method should be used, go to the section Selecting a Mounting Method.
- If some of the parameters have to be changed, go to the section Method Edit.


### Select a Mounting Method

The display should show the Method Group Menu:

```
SELECT METHOD GROUP ↓
1. STRUERS 2. GROUP NUMBER 2
3. GROUP NUMBER 3 4. GROUP NUMBER 4
F1:CONFIG. F3:MAN. FUNC.
```

If the heading in the display is different:

 Press Esc, until a display with the above heading appears.

↓  
 Press MENU UP/DOWN ▲ ▼ until you reach the desired Method Group, e.g. GROUP NUMBER 2.

↓  
 Press ENTER ↵, to accept the group.

↓

*ProntoPress-10*  
*Instruction Manual*

Group:GROUP NUMBER 2      SELECT METHOD ↓  
1. **METHOD 1**                      2.METHOD 2  
3. METHOD 3                              4.EMPTY METHOD  
F1: COPY                                  F3:RESET   F4:RENAME

↓  
 Press MENU UP/DOWN ▲ ▼ until you reach the desired Method, e.g. METHOD 3.

↓  
Group:GROUP NUMBER 2      SELECT METHOD ↓  
1. METHOD 1                              2.METHOD 2  
3. **METHOD 3**                              4.EMPTY METHOD  
F1: COPY                                  F3:RESET   F4:RENAME

↓  
 Press ENTER ↵, to see the method.

↓  
Grp: 2   Method: 3    Resin:**MultiFast**  
HEAT: Time: 7:00    T: 180°C   F: 20 kN  
COOLING:    Rate:    HIGH      Time: 3:00  
F1:SELECT PREHEATING


*Method Edit*

Select a Mounting Method. (The Struers methods in Group 1 are protected against editing).


```
Grp: 2 Method: 3 Resin:MultiFast
HEAT: Time: 7:00 T: 180°C F: 20 kN
COOLING: Rate: HIGH Time: 3:00
F1:SELECT PREHEATING
```

↓  
 Move the cursor, to the desired parameter to be changed, e.g. Heating Time.


↓  
Grp: 2 Method: 3 Resin:MultiFast  
HEAT: Time: **7:00** T: 180°C F: 20 kN  
COOLING: Rate: HIGH Time: 3:00  
F1:SELECT PREHEATING

↓  
 Press ENTER ↵, to allow editing of the parameter marked with [ ].

↓  
Grp: 2 Method: 3 Resin:MultiFast  
HEAT: Time: [7:00] T: 180°C F: 20 kN  
COOLING: Rate: HIGH Time: 3:00  
F1:SELECT PREHEATING

↓  
 Press MENU UP/DOWN   to increase/decrease the value.

↓  
Grp: 2 Method: 3 Resin:MultiFast  
HEAT: Time: [6:30] T: 180°C F: 20 kN  
COOLING: Rate: HIGH Time: 3:00  
F1:SELECT PREHEATING

↓  
 Press ENTER ↵ to accept the new value.

↓  
Grp: 2 Method: 3 Resin:MultiFast  
HEAT: Time: **6:30** T: 180°C F: 20 kN  
COOLING: Rate: HIGH Time: 3:00  
F1:SELECT PREHEATING F4:SAVE

If you have modified parameters in a method and switch off the voltage with this method open, the mounting press will start again in the same method, and with the modified parameters.

When modifications have been made, a SAVE option will appear as a function key option.

To save the modified method in the database: See Advanced Operations.

The default method is called: EMPTY METHOD. When modifications have been saved, the method is automatically renamed as: UNNAMED METHOD. You are recommended to edit the name for the methods you want to keep, for easier identification.

## **Placing the Specimen**

- Press and hold RAM UP ▲ to raise the lower ram to its upper limit.
- Apply "*Mould Release Agent*" to the surface of the lower ram.
- Place the specimen on the ram. The specimen must be clean, dry and free from grease. The distance between the specimen and the cylinder wall must be minimum 3 mm to avoid cracks in the resin.

## *Pouring Resin over the Specimen*

- Press and hold RAM DOWN ▼ to lower the ram to its lower limit.
- Fill a suitable amount of resin into the cylinder by means of the enclosed funnel.

**Always** make sure there is sufficient resin to cover the sample after compression. Please note that the volume of the resin reduces when the granulate becomes compressed. If insufficient resin is used, the rams may come in contact with the sample, and the cylinder and the rams may be damaged.

### ***Important***


The "*Mould Release Agent*" must always be applied to the mounting rams as a thin layer to prevent the resins sticking to the surface. By means of Struers AntiStick a thin layer of stearate powder can easily be dabbed on the rams.

### Placing Two Specimens

- Follow the instructions in Placing the Specimen.
- Apply "*Mould Release Agent*" to all surfaces of the intermediate ram.
- Place the intermediate ram on top of the resin.
- Place the second specimen on the intermediate ram and fill the cylinder with an adequate quantity of resin.

Again, always make sure there is sufficient resin to cover the sample after compression.

### Installing the Top Closure

- Remove resin dust from the upper part of the mounting cylinder.
- Clean the cylindrical surface of the upper ram. Cured resin can easily be removed without damage to the surface of the ram using the scraper supplied.
- Apply "*Mould Release Agent*" to all accessible surfaces of the upper ram.
- Place the top closure with the upper ram on the mounting cylinder.
- Press the top closure straight down, turning it counter-clockwise until you hear a click.
- Press the top closure down and turn it clockwise until its lower limit.
- Turn the top closure a quarter of a turn back.

***Important***

If the ram does not fit in the cylinder then check ram and cylinder for cured resin. The tolerance between the cylinder and the ram is very fine and even small amounts of resin from previous mounting may cause problems.

### **Starting the Mounting Process**

- Select a Mounting Method. See Select a Mounting Method.
- Press START  $\blacklozenge$ , and the process will run automatically.

### *Display during the Mounting Process*

The display with the mounting parameters will change to a new display, which will constantly monitor the mounting process:

MOUNTING PROCESS: HEATING	
Remaining Process Time:	8:30
Temperature (pre-set/act.):	180 / 145°C
Force (pre-set/act.):	20 / 20 kN

The time units are: minutes:seconds.

The parameters cannot be changed for the step being carried out.

### **Stopping the Mounting Process**

- The machine automatically stops when the cooling time has elapsed. The machine can be stopped at any time during the mounting process by pressing STOP  $\odot$ .

<p><i>If you have stopped the machine during the mounting process:</i></p> <p>Cool the mounting cylinder in minimum 2 min before opening, after a heating period. Please note that the mount possibly will be destroyed.</p>
--

### **Removing the Top Closure**

When it is indicated in the display, that the mounting process is finished:

- Turn the top closure counter-clockwise until released from thread.
- Press RAM UP  $\blacktriangle$  to raise the lower ram to its upper limit.
- Turn the top closure to one side to give access to remove the specimen.

# Reference Guide


Table of Contents	Page
<b>1. Advanced Operations</b>	
Installation.....	23
Recirculation Cooling Unit (Optional).....	23
Socket for Recirculation .....	23
Software Settings.....	24
Configuration Menu .....	24
Setting the Temperature Unit .....	25
Setting the Force or Pressure Mode .....	26
Setting the Force Unit.....	27
Setting the Pressure Unit .....	28
Setting the Countdown Mode .....	29
Mounting Method Options .....	30
Saving a Method.....	30
Copying a Method.....	31
Insert Method .....	32
Reset a Method .....	34
Editing Names .....	35
Name Editing Principles .....	35
Process Control .....	36
Preheating.....	36
Heating Time .....	37
Cooling .....	37
Manual Functions.....	38
Heating of Cylinder .....	38
Cooling of Cylinder .....	39

<b>2. Struers Metalog Guide™</b> .....	41
<b>3. Application Guide for Hot Mounting</b> .....	41
<b>4. Accessories</b> .....	42
<b>5. Consumables</b> .....	43
<b>6. Trouble-Shooting</b> .....	44
<b>7. Maintenance</b>	
Daily Service .....	48
After 200 mountings .....	48
Removing residues.....	48
Resetting the counter .....	48
Cleaning the Rams .....	48
Removing the Upper Ram .....	49
Lubricating the Threads of the Top Closure .....	49
Decalcifying the Cooling Coil .....	49
Replacing the Cooling Water.....	50
Checking the Recirculation Cooling Unit.....	50
<b>8. Technical Data</b> .....	51
<b>9. Menu Overview</b> .....	53

## 1. Advanced Operations

### Installation


#### Recirculation Cooling Unit (Optional)

A Recirculation Cooling Unit (TRECI) can be connected by means of a Recirculation Cooling Unit Connector (ROPRE).

Connect the cable from the pump in TRECI to ROPRE, according to the diagram of ROPRE.

#### Socket for Recirculation

- ② Outlet hose
- ⑤ Socket for recirculation

- Connect the cable with plug from ROPRE to the socket (pos. 5) on the back of the mounting press.
- Connect the cable from ROPRE to the mains.
- Check that the pump rotates in the direction indicated on the pump. If the direction is incorrect, switch two of the phases in the connection to the pump in TRECI (3-phase versions only).
- Replace the outlet tube on the pump in TRECI with the non-return valve with elbow and hexagon nipple from ROPRE. Struers recommends that you seal the threads.
- Connect the pressure hose of the mounting press to the hexagon nipple. Remember to mount the gasket.
- Lead the outlet hose (pos. 2) to the inlet on TRECI. Be absolutely sure to place the hose with a steady downward slope and without obstructions.
- Attach the water outlet to the inlet of TRECI.

**Software Settings**  
*Configuration Menu*

Switch on the power at the main switch. The following display will appear briefly:

```
ProntoPress-10  
Version: 2.0
```

The mounting press will start in the same state, as when the voltage was last switched off. The display to appear could be the Select Method Group Menu. If the heading in the display is different, press Esc until the Select Method Group Menu appears.

The Select Method Group Menu is the highest level in the menu structure. From here you can go to configuration, manual functions, select mounting methods in the database and open the mounting methods.

```
SELECT METHOD GROUP ↓  
1. STRUERS                2. GROUP NUMBER 2  
3. GROUP NUMBER 3          4. GROUP NUMBER 4  
F1:CONFIG.                F3:MAN. FUNC.
```

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




Press ENTER  to activate the Temperature Unit Menu.


Press MENU UP/DOWN  ▲ ▼ to select the Temperature Unit you prefer.

Press ENTER  to accept the Temperature Unit.

*Setting the Force or Pressure Mode*

FORCE mode is the common function mode for a mounting press, setting the actual force on the ram.  
In PRESSURE mode, you set the absolute pressure independently of the cylinder diameter.

 Press MENU UP/DOWN  to select the Force or Pressure Parameter.


↓  

CONFIGURATION		↓
Cylinder Diameter:	50 mm	
Temperature Unit:	CELSIUS	
Force or Pressure:	<b>FORCE</b>	Unit: kN

↓  
 Press ENTER  to activate the Force or Pressure Menu.

↓  

FORCE OR PRESSURE	
<b>FORCE</b>	
PRESSURE	

↓  
 Press MENU UP/DOWN  to select the mode you prefer.

↓  
 Press ENTER  to accept the mode.

*Setting the Force Unit*

The press has to be in Force mode. See above.

 Press MENU UP/DOWN  to select the Force Unit parameter.


↓  

CONFIGURATION		↓
Cylinder Diameter:	50 mm	
Temperature Unit:	CELSIUS	
Force or Pressure:	FORCE	Unit: <b>kN</b>

↓  
 Press ENTER  to activate the Force Unit Menu.

↓  

FORCE UNIT	
Lbf	
kp	
<b>kN</b>	

↓  
 Press MENU UP/DOWN  to select the unit you prefer.

↓  
 Press ENTER  to accept the Force Unit.

*Setting the Pressure Unit*

The press has to be in Pressure mode. See above.

 Press MENU UP/DOWN  to select the Pressure Unit Menu.


↓  

CONFIGURATION		↓
Cylinder Diameter:	50 mm	
Temperature Unit:	CELSIUS	
Force or Pressure:	PRESSURE	Unit: <b>MPa</b>

↓  
 Press ENTER  to activate the Pressure Unit Menu.

↓  

PRESSURE UNIT	
psi	
bar	
<b>MPa</b>	

↓  
 Press MENU UP/DOWN  to select the unit you prefer.

↓  
 Press ENTER  to accept the pressure unit.

*Setting the Countdown Mode*

The Countdown mode **COUNT FROM START** is the common control of heating time on a mounting press, setting the heating time, and allowing the mounting press to start and stop independently of the temperature.


In the Countdown mode **COUNT FROM PRE-SET TEMP.**, you set the countdown time for preheating and heating counting from the time, when the temperature is reached. In this case values for time should be reduced, compared to the countdown mode **COUNT FROM START**, and the heating process will be independent of the initial temperature of the mounting unit.

 Press MENU UP/DOWN  to select the Countdown Mode.

↓

 Press ENTER  to activate the Countdown Menu.

↓

 Press MENU UP/DOWN  to select the mode you prefer.

↓

 Press ENTER  to accept the Countdown mode.

Setting the Language and Cylinder Diameter: see the section Getting Started.
---

## Mounting Method Options

### *Saving a Method*

You can save the changes you have made, while working with a specific method, in the database.

- Edit the method. See Getting Started.

```
Grp: 2 Method: 3 Resin:MultiFast
HEAT: Time: 6:30 T: 180°C F: 20 kN
COOLING: Rate: HIGH Time: 3:00
F1:SELECT PREHEATING F4:SAVE
```

**F4** Press F4:SAVE.

```
Save changes in
Group: 2. GROUP NUMBER 2
Method: 3. METHOD 3
ENTER: YES ESC: NO
```

**↵** Press ENTER ↵ to save the modified method.

### ***Important***

When saving the changes, the original method will be overwritten. If you want to keep the original method, you should make a copy of the method with a new name, thus making the changes in a copy rather than changing the original method. See Copying a Method.

*Copying a Method*

Copying a method is a shortcut to creating a new mounting method on the basis of an existing one.

- Select the group with the mounting method you want to copy from. See Getting Started, Select a Mounting Method.


```
Group:GROUP NUMBER 2   SELECT METHOD ↓
1.METHOD 1           2.METHOD 2
3.METHOD 3         4.EMPTY METHOD
F1:COPY  F2:INSERT  F3:RESET
```

↓  
 Press ENTER ↵ if you like to see the method parameters.

↓  
 Press Esc to return to the above display.

↓  
 Press F1:COPY to copy the method to the buffer.

```
↓
Do you want to copy
Group:  2.  GROUP NUMBER 2
Method: 3.  METHOD 3      to buffer ?
ENTER: YES                      ESC: NO
```

↓  
 Press ENTER ↵ to accept copying the method to the buffer.


*Insert Method*

If you want to insert the method in a different group:

 Press Esc to return to the Method Group Screen.

↓  
 Press MENU UP/DOWN   to select the group, in which you want to insert the method.

↓  
 Press ENTER  to accept the group.

↓  
 Press MENU UP/DOWN   to select the method number, in which you want to insert the method.  
↓

If it is not an EMPTY METHOD:

Press ENTER ↵ to see the method.

Press Esc to return to the above display.

Press F2:INSERT to insert the method from the buffer.

```
Do you want to copy buffer to
Group:   3.  GROUP NUMBER 3
Method:  4.  EMPTY METHOD      ?
ENTER: YES                               ESC: NO
```

Press ENTER ↵ to accept inserting the method.

**Important**

When you insert to a new method, the existing parameters in this method will be overwritten.  
However, if the method has not been used, or if it has been reset, the method name will be: EMPTY METHOD.

*Reset a Method*

If a method is not relevant any more, it should be reset.

The parameters will change to default values, which can easily be changed to a new method later on.

The method name will change to: **EMPTY METHOD**, telling you, that you can copy to or modify in this method without replacing any valuable data.

- Select the group with the mounting method you want to reset.
- Select the mounting method you want to reset.

```
Group:GROUP NUMBER 2   SELECT METHOD ↓
1.METHOD 1           2.METHOD 2
3.METHOD 3           4.EMPTY METHOD
F1: COPY   F2: INSERT  F3: RESET
```

If it is not an **EMPTY METHOD**:

Press ENTER ↵ to see the method.

Press Esc to return to the above display.

Press F3:RESET.

```
Do you want to reset this method ?
Group:   2. GROUP NUMBER 2
Method:  3. METHOD 3
ENTER: YES                               ESC: NO
```

Press ENTER ↵ to reset the method.

*Editing Names*

The names for the method groups, methods and mounting resins can be edited and changed to suit your preference.



```
Group:GROUP NUMBER 2   SELECT METHOD ↓
1.METHOD 1           2.METHOD 2
3.METHOD 3           4.EMPTY METHOD
F1:COPY   F2:INSERT F3:RESET F4:RENAME
```

↓  
**F4** Press F4:RENAME.


↓  
Text:EMPTY METHOD  
ABCDEFGHIJKLMNOPQRSTUVWXYZ 0123456789  
, . ; : - ! " # % & / ( ) = ? +  
F1:LEFT F2:SPACE F3:RIGHT F4:DELETE

*Name Editing Principles*

- Place the main cursor on the character you want to change, using F1:LEFT or F3:RIGHT. Use MENU DOWN  to move to the character set in the next line. An auxiliary cursor in the text line shows the position in the method name.
- Write the new name using the following keys:
  - F1 Moves the main cursor to the left
  - F3 Moves the main cursor to the right
  - F2 Inserts a space in the text
  - F4 Deletes one character in the text
  - ↵ ENTER places the new character in the method name and moves the auxiliary cursor in the name to the right. Repeat the procedure for each character.

 Press MENU DOWN  to move to the character set.

↓ Write the new name using the above name editing procedures.

 Press MENU UP  to return to the text line with the method name.

↓  
 Press ENTER ↵ to accept the new name.

## Process Control

### *Preheating*

A preheating phase can improve the quality of the mount in some applications.

- To activate the preheating in a mounting method:

```
Grp: 2 Method: 3 Resin:MultiFast
HEAT: Time: 7:00 T: 180°C F: 20 kN
COOLING: Rate: HIGH Time: 3:00
F1:SELECT PREHEATING
```

**F1** Press F1:SELECT PREHEATING to activate preheating.

```
Grp: 2 Method: 3 Resin:MultiFast
PRE : Time: 0:30 T: 180°C F: 0 kN
HEAT: Time: 7:00 T: 180°C F: 20 kN
F1:BYPASS PREHEATING
```

- To deactivate the preheating:

**F1** Press F1:BYPASS PREHEATING to deactivate preheating.

### Heating Time


The machine can count the preheating and heating time either from start or from when the pre-set temperatures are reached.

To set the Counting mode for heating time: See Software Settings, Setting the Countdown Mode.

### Cooling

The cooling time is always counted from start of cooling.

The machine can count through the cooling with all three different cooling rates.

**ECONOMY** cooling rate will only allow 20 % of the water flow during the whole cooling period, compared to fully open circulation (cooling rate: HIGH). The cooling time will be slightly longer than for cooling rate: HIGH.

**LOW** cooling rate will only allow 2% of the water flow in the first 3 minutes of the cooling period, then 20% of full flow. This will allow for a very uniform cooling of the mount, thus reducing the internal stress in the resin. The cooling time will in some cases be slightly longer than for cooling rate: HIGH.

**HIGH** cooling rate will give full flow during the whole cooling period. This gives a very rapid cooling.

**Manual Functions**  
*Heating of Cylinder*

The mounting cylinder can be heated without pressure on. It can be closed or open. This feature can be used as stand by, to reduce the process time, or for special applications.

Go to the Select Method Group Menu:

```
SELECT METHOD GROUP ↓
1. STRUERS          2. GROUP NUMBER 2
3. GROUP NUMBER 3    4. GROUP NUMBER 4
F1:CONFIG.          F3:MAN. FUNC.    F4:RENAME
```

If the heading in the display is different, press Esc until the Select Method Group Menu appears.

**F3** | Press F3:MAN. FUNC to activate the manual functions.

```
↓
State: OFF          HEATING OF CYLINDER ↓
Pre-set Temperature:    50 °C
Actual Temperature:    25 °C
F1:START HEATING
```

Adjust the pre-set temperature to the desired value. See Method Edit.

For safety reasons, the temperature can only be selected in the range 30 to 60°C. The heating is automatically switched off after 1 hour.

To start heating:

**F1** | Press F1:START HEATING and state changes to ON.

To stop heating:

**F1** | Press F1:STOP HEATING and State changes to OFF.

The heating is automatically switched OFF, if manual cooling is selected or if a mounting process is started.

### *Cooling of Cylinder*

The mounting cylinder can be cooled without pressure on. This can be performed with the top closure removed or in place. This feature can be used for supplementary cooling, if, for example, an ongoing mounting process should be stopped.

Go to the Select Method Group Menu:

```
SELECT METHOD GROUP ↓
1. STRUERS           2. GROUP NUMBER 2
3. GROUP NUMBER 3     4. GROUP NUMBER 4
F1:CONFIG.           F3:MAN. FUNC.  F4:RENAME
```

If the heading in the display is different, press Esc until the Select Method Group Menu appears.

**F3** Press F3 to activate the manual functions.

↓

```
State: OFF           HEATING OF CYLINDER ↓
Pre-set Temperature: 50 °C
Actual Temperature:  25 °C
F1:START HEATING
```

↓

 Press MENU UP/DOWN  to come to the display for cooling:

↓

```
State: OFF           COOLING OF CYLINDER ↑
Cooling Time:       4:30
Cooling Rate:       HIGH
F1:START COOLING
```

- Adjust the cooling time and cooling rate to the desired values. See Method Edit.
- If the cooling time is set to 0:00, the cooling time will automatically count from the maximum value: 15:00, when the cooling is started.

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The cooling is automatically switched off, if manual heating is selected or if a mounting process is started.

## **2. Struers Metalog Guide™**

Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method, including choice of consumables is easy. Always consult Struers Metalog Guide™ for the correct preparation method for the actual specimens.

Struers Metalog Guide™ contains 6 useful chapters:

*Metalogram*

A quick and safe guide to the right preparation method.

*Metalog Method*

A complete catalogue of preparation methods, based on Struers' extensive materialographic experience, and employing Struers' range of consumables.

*Preparation Philosophy*

The basics of modern specimen preparation, seen from a professional point of view.

*Metalog Process*

The materialographic preparation process from start to finish, logically explained.

*Metalog Master*

A combined trouble-shooting guide and supply of in-depth information on the processes of mechanical preparation, including an expert system for the solving of preparation problems.

*Metalog Code*

Quick access to the relevant consumables for the chosen preparation methods.

### ***Struers Metalog Guide™***

A complete guide to materialographic specimen preparation. Contact your local dealer for a free copy of Metalog Guide™.

## **3. Application Guide for Hot Mounting**

Useful mounting data and hints can be found in the Struers [Application Guide for Hot Mounting](#) or visit the Struers website on <http://www.struers.com> and see the section on **Knowledge**.

## 4. Accessories

<b>Type</b>	<b>Description</b>	<b>Code</b>
<b>Mounting Unit</b> consisting of cylinder, upper and lower ram, heating/cooling unit and top closure	25 mm dia. 30 mm dia. 1 ¼" dia. 1 ½" dia. 40 mm dia. 50 mm dia.	CYLM1 CYLTE CYLQU CYLHA CYLIF CYLAM
<b>Intermediate Ram</b> for making 2 simultaneous mounts	25 mm dia. 30 mm dia. 1 ¼" dia. 1 ½" dia. 40 mm dia. 50 mm dia.	CYLDO CYLDO CYLDO CYLDO CYLDO CYLDO
<b>Recirculation Cooling Unit</b>		TRECI
<b>Recirculation Cooling Unit Connector</b> for connection of Recirculation Cooling Unit (TRECI) to grinding/polishing machines		ROPRE

## 5. Consumables

Resin	Application	Specific property	Material	Contents	Code
<b>ConduFast</b> (Resin 1)	Electropolishing	Electrically conductive Very low shrinkage	Acrylic resin with iron powder filler Thermoplastic, Grey	1 kg	RESFE
<b>SpeciFast</b> (Resin 3)	Glass clear mounts Porous samples Surface electrical insulator for ConduFast	Transparent Medium shrinkage	Acrylic Thermoplastic	1.5 kg 7.5 kg 25 kg	RESTH RESYV RESFA
<b>PolyFast</b>	Edge retention	Very low shrinkage Medium removal rate	Phenolic resin with carbon filler Thermosetting, Black	1 kg 7.5 kg	FAPSA FAPME
	Examination in scanning electron microscope	Electrically conductive Low emission in the vacuum chamber.			
<b>IsoFast</b> (Resin 4)	Edge retention Planeness	Low shrinkage Low removal rate	Dialylphthalate with glass fibre filler Thermosetting, Green	1 kg 7.5 kg	RESDI RETTY
<b>DuroFast</b> (Resin 5)	Edge retention Planeness	Very low shrinkage Good adhesion Very low removal rate	Epoxy with mineral filler Thermosetting, Black	1 kg 7.5 kg	RESIF RESEN
<b>MultiFast</b>	Routine examination, Backup resin	Low shrinkage Medium removal rate	Phenolic resin with wood flour filler Thermosetting, Black	2.5 kg 7.5 kg 25 kg	RESRU RESUK RESLA
<b>MultiFast Green</b>	Routine examination Colour marking	Low shrinkage Medium removal rate	Phenolic resin with wood flour filler Thermosetting, Green	2.5 kg 7.5 kg 25 kg	FAGSA FAGME FAGLA
<b>MultiFast Brown</b>	Routine examination Colour marking	Low shrinkage Medium removal rate	Phenolic resin with wood flour filler Thermosetting, Brown	2.5 kg 7.5 kg 25 kg	FABSA FABME FABLA

Product	Application	Specific property	Material	Contents	Code
<b>Pre-Mount</b>	Serial mounting of uncomplicated shapes	Easy to handle Low shrinkage Medium removal rate	Phenolic resin with wood flour filler In tablet form Thermosetting, Black	<b>25 mm/ 1" dia.</b> approx. 1100 pcs. 10 kg	RESUN
				<b>30 mm/ 1.1/4" dia.</b> approx. 760 pcs. 10 kg	RESGI
				<b>40 mm/ 1 1/2" dia.</b> approx. 460 pcs. 10 kg	RESAM

Product	Application	Material	Contents	Code
<b>Mount Release Wax</b>	To avoid sticking of the mounts to the rams	Wax	24 g	PREWA
<b>AntiStick</b>	To avoid sticking of the mounts to the rams	Very fine stearate powder	8 g	FASTI

## 6. Trouble-Shooting

Display/error	Cause	Action
<b>Error Message</b>		
Mounting unit not connected	The plug for the mounting unit is not connected properly. (Appears when START $\diamond$ is pressed).	Insert the plug firmly. If message still appears, call a Struers service technician.
PT100 amplifier uncalibrated	The temperature measurement mechanism is uncalibrated or there is an error in calibration. Appears during start-up mode together with a long beep.	Call a Struers service technician.
Pressure Error	Defect in the pressure system. Appears, when the pre-set pressure is not reached within 25 seconds.	Call a Struers service technician.
Mounting pressure too high	Formation of gas in the mounting cylinder. Appears, if the actual force should exceed 65 kN.	Cool the mounting unit down completely. Move the lower ram down. Open the top closure.
Error, hydraulic-pump overloaded	The motor has been overloaded. Appears, either when pressing: RAM UP/DOWN $\blacktriangle$ $\blacktriangledown$ or during the process.	Let the motor cool down and start the mounting press again. If it does not help: Call a Struers service technician.
Heating Error	The plug for the mounting unit is not correctly connected. Appears, when the pre-set temperature is not reached within 6 minutes.	Re-connect the plug.
	Defect in the heating system.	Call a Struers service technician.
Insufficient Cooling	Mains tap supplying cooling water is either closed or not opened sufficiently. Filter at the water inlet is blocked. The water is too hot. Appears, when the cooling is inadequate.	Open the tap. Clean the filter.
	Insufficient water within the Recirculation Cooling Unit. The water is too hot. Appears, when the cooling is inadequate.	Fill to the correct water level. See the section: Maintenance.
	Defect in the cooling system.	Call a Struers service technician.
Top closure not unscrewed	The top closure has not been unscrewed completely after a mounting. Appears when RAM UP $\blacktriangle$ is pressed after a mounting.	Turn the top closure until it is completely free of the threads. If this has already been done, the problem is that the closure cannot spring free. Hold the swivel arm down with one finger and press RAM UP $\blacktriangle$ again.

*ProntoPress-10*  
*Instruction Manual*

<b>Display/error</b>	<b>Cause</b>	<b>Action</b>
<b>Warning</b>		
Top closure not tightened	Top closure not mounted correctly. Appears, when pressing: START ◇.	Screw the top closure to the bottom position and unscrew it ¼ turn. If this does not help: Call a Struers service technician.
Pre-set force too high	The pre-set value of force exceeds the maximum force on the piston rod. This can happen if the configuration has been changed to a larger cylinder diameter.	No action is needed <i>Alternatively:</i> if you want to bring the force within the maximum force of the press: decrease the value.
Pre-set pressure too high	The pre-set value of pressure exceeds the maximum force on the piston rod. This can happen if the configuration has been changed to a larger cylinder diameter.	No action is needed <i>Alternatively:</i> if you want to bring the pressure within the maximum force of the press: decrease the value.
<b>Acoustic Signals</b>		
Long beep.	The command can not be accepted.	Check User's Guide, Section 2.
Four double beeps.	There is an error.	See the error message.

*ProntoPress-10*  
*Instruction Manual*

<b>Display/error</b>	<b>Cause</b>	<b>Action</b>
<b>Machine Problems</b>		
Insufficient compression.	Incorrect setting of force/pressure.	Set the correct parameter. The force should be min. 5 kN/1000 Lbf/500kp.
	Incorrect configuration of cylinder diameter.	Check the configuration.
	Incorrect unit for force or pressure.	
Insufficient heating.	Incorrect setting of time for preheating or heating.	Set the correct parameter.
	Incorrect temperature unit set.	Check the configuration.
Insufficient cooling.	Incorrect setting for cooling time.	Set the correct parameter.
	Incorrect setting for cooling rate.	
	Incorrect Temperature unit set.	Check the configuration.
	Mains tap supplying cooling water is either closed or not opened sufficiently. Filter at the water inlet is blocked. The water is too hot.	Open the tap. Clean the filter.
	Insufficient water within the Recirculation Cooling Unit. The water is too hot.	Fill to the correct water level. See the section: Maintenance.
	Defect in the cooling system.	Call a Struers service technician.
Cooling water drips underneath the machine.	The quick coupling is not mounted correctly.	Remove the cover from the mounting unit and check the quick coupling connections.
The top closure will not start threading on the mounting cylinder.	The top closure is not mounted correctly.	Press the top closure straight down, turning it counter-clockwise until you hear a click. Turn the top closure clockwise.
	The upper ram is too hot.	Let the top closure and the upper ram cool down. Reduce the mounting temperature.
	Cured resin in the upper part of mounting cylinder.	Clean the mounting cylinder with a brass wire brush.
	Cured resin on the cylindrical surface of the upper ram.	Clean the ram with the scraper supplied.
	The top closure has been dropped on the floor/table, causing a bulge in the edge of the ram.	Replace the upper ram.
	The swivel arm has become askew.	Call a Struers service technician. In the meantime the mounting press can be used without the swivel arm.

*ProntoPress-10*  
*Instruction Manual*

Display/error	Cause	Action
	Damage to the threaded connection in the top closure or to the upper ram.	Remove the upper ram from the top closure (see instructions in chapter Maintenance). Try both the following procedures to identify the problem: <ul style="list-style-type: none"> <li>■ Try to mount the top closure, without the upper ram. If this is not possible, call a Struers service technician.</li> <li>■ Try to insert the disconnected upper ram into the mounting cylinder. If this is not possible, call a Struers service technician.</li> </ul>
The top closure cannot be screwed down completely.	Dirt in the threads on the top closure and the mounting cylinder.	Clean the threads. Use only a dry lubricant powder.
	The disc for thermal insulation, placed on the top of the upper ram, has a larger diameter than the upper ram.	Call a Struers service technician.
The top closure can not be loosened.	<p>The top closure has not been turned ¼ turn back before starting the process.</p> <p>Cured resin on the cylindrical surface of the top ram.</p> <p>Dirt in the threads in the top closure</p>	<p>To release the top closure: Move the lower ram up and down several times.</p> <p>If this does not help: -Put heating on for 1 min.</p> <p>If this does not help: - Set the force or pressure to zero. - Set the heating time and cooling time to 15 min. - Complete a mounting process.</p> <p>If this does not help: - Remove the two handles on the top closure. - Remove the plastic cover from the top closure (See Maintenance). - Loosen the top closure with a fork spanner.</p>
JAPANESE TEXTS	The display can show the screen texts in Japanese. When choosing a European language in the LANGUAGE menu, you might accidentally press ENTER ↵ at the Japanese option.	To return to English or French texts, just follow this procedure: <ul style="list-style-type: none"> <li>■ Press Esc until a long beep is heard.</li> <li>■ Press F1.</li> <li>■ Press MENU  four times.</li> <li>■ Press ENTER ↵.</li> </ul>

## 7. Maintenance

### Daily Service

- Clean all accessible surfaces with a moist cloth.

### After 200 mountings

When mounting, some old cured resin will fall down from the lower ram and built up underneath. The cured resin must be cleaned away, to avoid it harming the machine. The axle journal securing the lower ram has been designed to break in this case, thus preventing serious damage to the machine.

After 200 mountings a message will appear in the display, reminding you to remove the resin residues.


### *Removing residues*

- Remove the mounting unit. See this section.
- Remove the resin under the lower ram with a cloth or a soft brush.
- Install the mounting unit. See this section.

### *Resetting the counter*

- Remove the message and reset the counter by pressing F1 when the message appears.

### Cleaning the Rams

Before each mounting, check that there is no resin left on the flat surface of the rams from previous operations. Clean the cylindrical surface of the upper ram and the intermediate ram (if fitted). Cured resin can easily be removed without any damage to the surface of the rams using the scraper supplied.


If the rams have been scratched on the polished surfaces, disassemble and polish them with 6 or 3  $\mu\text{m}$  diamonds on a polishing cloth.

#### **Important**

Do not polish the surface on the sides of the rams. If the sides of the rams are heavily scratched, they must be scrapped. If the top closure has been dropped, causing a dent or deformation in the edge of the upper ram, the ram must also be scrapped.

## Removing the Upper Ram

- Take the cover plugs out of the ends of the handles.
- Unscrew the two handles with an Allen key (pos. 1).
- Remove the plastic cover (pos.2).
- Lift the upper ram and remove the lock washer (pos.3).
- Pull out the upper ram. Do not remove the rod (pos. 4) from the upper ram unless absolutely necessary.
- If it is necessary to remove the rod, hold the upper ram in a vice or a similar device. The ram surface **must** be protected with plastic or soft metal.

- ① Handle
- ② Plastic cover
- ③ Lock washer
- ④ Rod

## Lubricating the Threads of the Top Closure

In case of a slight friction in the threads, the lubrication should be carried out using a dry lubricant such as MoS<sub>2</sub> Lubricant Powder or graphite.

### **IMPORTANT**

Please observe that the threads in the joint between the top closure and the mounting cylinder must not be lubricated with oil or grease. As the mounting cylinder operates at high temperatures, only a dry lubricant, resistant to high temperatures, can be used.

## Decalcifying the Cooling Coil

When using cooling water from the mains tap in areas with a high chalk or mineral content, deposits can build up in the cooling coil. This may reduce the cooling effect, so once a year the cooling coil should be decalcified.

- Remove the mounting unit. (See this section.)
- Flush the coil with a **mild** decalcifying fluid, such as that used for coffee machines.
- Flush the cooling coil with clean water.
- Re-install the mounting unit.

## Replacing the Cooling Water

*Only relevant when a Recirculation Cooling Unit is connected:*


- Replace the cooling water in the Recirculation Cooling Unit at least once a month.
- Place a Disposable Plastic Container (TREPO) in the tank and fold it over the edge.
- Fill the tank with 30 l of water. NB! Too high a level of coolant in the tank might damage the pump.
- Add Struers Additive according to the instructions.
- Put the lid on the tank.
- Place the pump in the lid

### **IMPORTANT**

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

## Checking the Recirculation Cooling Unit

The cooling unit should be checked for cooling water every week and refilled when the flushing pump cannot reach the cooling water. Remember to use Struers Additive (pH should be around 9). Replace if pH is  $\leq 8.0$ . To check the concentration of Additive use a refractometer (0-10 Brix, see diagram.). Concentration = 2 x Brix.

## 8. Technical Data


Subject		Specifications	
		Metric/International	US
Mounting Units (Optional)	Diameter	25, 30, 40, 50 mm	1 ¼", 1½"
Compression	Force on piston rod	<i>Preheating:</i> 0-20 kN in steps of 5 kN 0-2000 kp in steps of 500 kp <i>Heating:</i> 0 - 50 kN in steps of 5 kN 0 - 5000 kp in steps of 500 kp	<i>Preheating:</i> 0 - 5000 Lbf in steps of 1000 Lbf <i>Heating:</i> 0 - 12000 Lbf in steps of 1000 Lbf
	Pressure in mounting cylinder	The pressure range is dependent on the cylinder diameter: <i>Preheating:</i> 0-10 to 40 MPa in steps of 2.5 to 10 MPa 0-100 to 400 bar in steps of 25 to 100 bar <i>Heating:</i> 0-25 to 100 MPa in steps of 2.5 to 10 MPa 0-250 to 1000 bar in steps of 25 to 100 bar	The pressure range is dependent on the cylinder diameter: <i>Preheating:</i> 0-1500 to 6000 psi in steps of 250 to 1000 psi <i>Heating:</i> 0-3750 to 15000 psi in steps of 250 to 1000 psi
Preheating	Temperature	100 - 180 °C in steps of 10 °C	212 - 356 °F in steps of 18 °F
	Time	0:30 - 15:00 min in steps of 0:30 min	
Heating	Temperature	<i>Process:</i> 100 - 180 °C in steps of 10 °C <i>Manual heating:</i> 30 - 60 °C in steps of 10 °C	<i>Process:</i> 212 - 356 °F in steps of 18 °F <i>Manual heating:</i> 86 - 140 °F in steps of 18 °F
	Time	0:30 - 15:00 min in steps of 0:30 min	
Cooling	Rate	<i>Economy</i> Intermittent flow: 20 % of full flow <i>Low</i> Intermittent flow: 2% of full flow for first 3 minutes Afterwards: 20% of full flow <i>High</i> Full flow	
	Time	0 - 15:00 min in steps of 0:30 min	

*ProntoPress-10*  
*Instruction Manual*

Subject		Specifications	
		Metric/International	US
Software and Electronics	LC display	4 x 40 characters	
	Controls	Touch pad	
	Memory	EPROM, RAM, and non-volatile RAM 8 kb containing 100 mounting methods	
	Socket for connected equipment	For connector (ROPRE) to be connected with Recirculation Cooling Unit (TREC1)	
Environment	Safety Standards	General 89/392/EEC & 91/386/EEC & 93/44/EEC EN292-1, EN292-2, EN60204-1 (IEC 204-1) EMC 89/336/EEC & 92/231/EEC EN50081-2, EN50082-2 FCC part 15, subpart B	
	Noise level	Approx. 60 dB(A), at a distance of 1.0 m/39.4" from the front of the machine	
	Surrounding temperature	5-40°C/41-104°F	
	Humidity	0-95% RH non condensing	
Water Supply	Tap Water or Recirculation Cooling Unit		
	Pressure for tap water	1 - 10 bar	14.5 - 145 psi
	Inlet	1/2" or 3/4"	1/2" or 3/4"
Supply Voltage	Power phases	1	
	Power consumption: Max. Total Mounting Unit	1500W 750W	
	Voltage/frequency:	Max. continuous load	
	100-120V / 50Hz	10A (100V) - 12A (120V)	
	100-120V / 60Hz	10A (100V) - 12A (120V)	
220-240V / 50Hz	4.5A (220V)- 5.0A (240V)		
220-240V / 60Hz	4.5A (220V)- 5.0A (240V)		
Dimensions and Weight	Width	420 mm	
	Depth	580 mm	
	Height (Mounting Unit and Top Closure installed)	420 mm	
	Weight	38 kg	

## 9. Menu Overview

# Quick Reference Card

## Placing the Specimen

- Press and hold RAM UP ▲ to raise the lower ram to its upper limit.
- Apply "Mould Release Agent" to the surface of the lower ram.
- Place the specimen on the ram.

## Pouring Resin over the Specimen

- Press and hold RAM DOWN ▼ to lower the ram to its lower limit.
- Fill a suitable amount of resin into the cylinder through the funnel.

## Installing the Top Closure

- Remove resin dust from the upper part of the mounting cylinder.
- Clean the cylindrical surface of the upper ram. Cured resin can easily be removed without damage to the surface of the ram using the scraper supplied.
- Apply "Mould Release Agent" to all accessible surfaces of the upper ram.
- Place the top closure with the upper ram on the mounting cylinder.
- Press the top closure straight down, turning it counter-clockwise until you hear a click.
- Press the top closure down and turn it clockwise until its lower limit.
- Turn the top closure a quarter of a turn back.

## Starting the Mounting Process

- Press START ⬠, and the process will run automatically.

## Stopping the Mounting Process

- The machine automatically stops when the cooling time has elapsed.

## Removing the Top Closure

- Turn the top closure counter-clockwise until released from the threads.
- Press RAM UP ▲ to raise the lower ram to its upper limit.
- Turn the top closure to one side to give access to remove the specimen.

## Select a Mounting Method

- Press Esc, until display with Method Group appears.
- Press MENU UP/DOWN ▲ ▼ until the desired Method Group appears.
- Press ENTER ↵, to accept the group.
- Press MENU UP/DOWN ▲ ▼ until the desired Method appears.
- Press ENTER ↵, again to see/edit the method.

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Date of Release 10.05.2005

## Spare Parts and Diagrams

*ProntoPress-10/-20*  
*Spare Parts and Diagrams*

---

**Always state *Serial No* and *Voltage/frequency*  
if you have technical questions or when ordering spare parts.**

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### **Table of contents**

### **Drawing**

#### **Common Parts to ProntoPress-10 and ProntoPress-20**

##### *Drawings*

ProntoPress-10, complete .....	15060002B
Bottom, Complete, ProntoPress-10.....	15060051Q
Cylinder Unit .....	15360036C
Top Closure .....	15360026A

#### **ProntoPress-10**

##### *Diagrams*

Hydraulic Diagram .....	15061002D
Wiring Diagram .....	15063450H

#### **ProntoPress-20**

##### *Diagrams*

Hydraulic Diagram .....	15071003C
Wiring Diagram .....	15073450F

***Please note***

Some of the drawings may contain **position numbers**  
not used in connection with this manual.

*ProntoPress-10/-20*  
*Spare Parts and Diagrams*

Drawing	Pos.	Spare Part	Cat. No.
15060002		<b>ProntoPress-10, complete</b>	
	50	Top cabinet, Left	15060112
		Top cabinet, Right (ProntoPress-20)	15070112
15060051		<b>Bottom, complete</b>	
	30+40	Rubber foot with screw and spacer	R4910140
	145	Fuse 3.15A T, glass tube	359MP109A
	350	Axle journal	15362907
	450	Coupling 45°	RNH11445
	510	Straight coupling with hose nipple	2NH10014
		Air Filter (Lid for hydraulic oil)	2YF90005
		<b>Various parts</b>	
		Funnel ø80/15	377MP096
		Measuring spoon	50300000
		Pipe Wrench 13mm	2GR00006
		Allen Key 5 (5 pcs)	RGR00050
		Pressure water hose	359MP001
	Gasket with filter 3/4"	445MP068	
	Gasket ø11/ø24x1.5	R3590359	
	Reducing piece 1/2" 3/4"	2NG30013	
15360036		<b>Cylinder Unit ø25</b>	
	155	Cover plate ø25	15368005
	190	Fixation screw	15050040
		Lower ram, ø25	R5368010
		Piston pin, for ø25 ram	R5360185
		Dust protection ring	15360313
		Scraper (PROAN)	05366901



*ProntoPress-10/-20*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Pos.</b>	<b>Spare Part</b>	<b>Cat No.</b>
15360026		<b>Top closure ø25</b>	
		Top closure ø25, complete	15360020
	30	Insulation, disc for ø25 top closure	15360300
	40	Pressure spring ø12.5x1.0	2GF10125
	50	Top part of ø25 top closure	15360170
	60	Lock washer	15050241
	70	Insulation shell for ø25 top closure	15360180
	80	Handle for top closure	15360027
15360036		<b>Cylinder unit ø30</b>	
	155	Cover plate ø30/1.1/4"	15368006
		Lower ram ø30	R5368011
		Piston pin, for ø30/ø1.1/4" ram	R5360186
15360026		<b>Top Closure ø30</b>	
		Top Closure ø30, complete	15360021
	30	Insulation disc for ø30/ø1.1/4" top closure	15360301
	50	Top part of ø30/1.1/4" top closure	15360171
	70	Insulation shell for ø30/ø1.1/4 top closure	15360181
15360036		<b>Cylinder unit ø1.1/4"</b>	
		Lower ram, ø1.1/4"	R5368012
15360026		<b>Top closure ø1.1/4"</b>	
		Top closure ø1.1/4", complete	15360022
15360036		<b>Cylinder unit ø1.1/2"</b>	
		Lower ram ø1.1/2"	R5368013
	155	Cover plate ø40/1.1/2"	15368007
		Piston pin, for ø40/ø1.1/2" ram	R5360187
15360026		<b>Top closure ø1.1/2"</b>	
		Top closure ø1.1/2", complete	15360023
	30	Insulation disc for ø40/ø1.1/2"/ø50 top closure	15360302
	50	Top part of ø40/1.1/2" top closure	15360172
	70	Insulation shell for ø40/ø1.1/2 top closure	15360182


*ProntoPress-10/-20*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Pos.</b>	<b>Spare Part</b>	<b>Cat No.</b>
15360036		<b>Cylinder unit ø40</b> Lower ram ø40	R5368014
15360026		<b>Top closure ø40</b> Top closure ø40, complete	15360024
15360036	155	<b>Cylinder unit ø50</b> Cover plate ø50 Lower ram ø50 Piston pin, for ø50 ram	15368008 R5368015 R5360188
15360026		<b>Top closure ø50</b> Top closure ø50, complete	15360025
	50	Top part of ø50 top closure	15360173
	70	Insulation shell for ø50 top closure	15360183

	Matr.:	Overfl.beh.: Nej	Målforhold:	Ikke ang. tol. efter DS/ISO 2768-		
		Projektionsmetode	1:5		Date	Sign.
B: Pos.80: MC5x10 fjernet. 301197 BMJ/					Tegn.	260594
A: Pos.80+90 tilf. (2-Uls -> 2-MC) 170596 BMJ/JTV	ProntoPress-10, komplet		Erst.:	15060002B		



- 450 510 Smøres med Unisilikon TK M 1011
- 60 340 445 Omnifit 100M
- 80 90 Limes med 409 GEL
- 90 High vacuum grease på 2 udv. flader.
- 150 Opsides. Flap placeres inde i hydraulik-rum.
- 480 Placeres 360mm fra kobling.

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




1. Hydraulic Powerpack
2. Double-acting Hydraulic Cylinder
3. ø5 Polyamid Tube
4. Pressure Transmitter
5. Valve Manifold
6. Relief Valve
7. Accumulator
8. Throttle Plug

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
SEE CONNECTIONS  
 DIAGRAM 15083150  
 WIRING FOR 120V AND 240V


- X1:1 RD2 9.5V IN1
- X1:2 RD2 9.5V IN2
- X1:3 RD3 24V IN1
- X1:4 RD4 24V IN2

- X7:1 NC.
- X7:2 NC.
- X7:3 W4 BK6 OUT ON
- X7:4 W4 BK5 N220V HEATING ON
- X7:5 W5 BK2 OUT ON
- X7:6 W5 BK1 N220V PUMP MOTOR
- X7:7 W3 BK4 M220V POWER IN
- X7:8 W3 BK3 N220V

- X6:1 W17 BL N24V SIGNAL FOR RECIRCUL. ON
- X6:2 W11 RD1 N24V RAM UP Y1
- X6:3 W12 RD1 N24V RAM DOWN Y2
- X6:4 W13 RD1 N24V COOLING Y3
- X6:5 W17 BK OUT1
- X6:6 W11 RD2 OUT2
- X6:7 W12 RD2 OUT3
- X6:8 W13 RD2 OUT4
- X3:1 W4 RD7 R IN THERMISTOR PT100
- X3:2 W14 RD1 +12V PRESSURE TRANSMITTER (BN) \*
- X3:3 W14 RD2 I IN TRANSMITTER (YE) (WH)
- X3:6 W14 AVREF
- X3:7 N.C.
- X3:8 RD5 +5V TOP CLOSED
- RD6 I IN NAMUR

\* BP1; ONLY VOLTAGE OUT TRANSDUCER.

REV. 1: BK - BK REV. 2: BK - BK REV. 3: BK - BK REV. 4: BK - BK REV. 5: BK - BK REV. 6: BK - BK REV. 7: BK - BK REV. 8: BK - BK REV. 9: BK - BK REV. 10: BK - BK REV. 11: BK - BK REV. 12: BK - BK REV. 13: BK - BK REV. 14: BK - BK REV. 15: BK - BK REV. 16: BK - BK REV. 17: BK - BK REV. 18: BK - BK REV. 19: BK - BK REV. 20: BK - BK REV. 21: BK - BK REV. 22: BK - BK REV. 23: BK - BK REV. 24: BK - BK REV. 25: BK - BK REV. 26: BK - BK REV. 27: BK - BK REV. 28: BK - BK REV. 29: BK - BK REV. 30: BK - BK REV. 31: BK - BK REV. 32: BK - BK REV. 33: BK - BK REV. 34: BK - BK REV. 35: BK - BK REV. 36: BK - BK REV. 37: BK - BK REV. 38: BK - BK REV. 39: BK - BK REV. 40: BK - BK REV. 41: BK - BK REV. 42: BK - BK REV. 43: BK - BK REV. 44: BK - BK REV. 45: BK - BK REV. 46: BK - BK REV. 47: BK - BK REV. 48: BK - BK REV. 49: BK - BK REV. 50: BK - BK REV. 51: BK - BK REV. 52: BK - BK REV. 53: BK - BK REV. 54: BK - BK REV. 55: BK - BK REV. 56: BK - BK REV. 57: BK - BK REV. 58: BK - BK REV. 59: BK - BK REV. 60: BK - BK REV. 61: BK - BK REV. 62: BK - BK REV. 63: BK - BK REV. 64: BK - BK REV. 65: BK - BK REV. 66: BK - BK REV. 67: BK - BK REV. 68: BK - BK REV. 69: BK - BK REV. 70: BK - BK REV. 71: BK - BK REV. 72: BK - BK REV. 73: BK - BK REV. 74: BK - BK REV. 75: BK - BK REV. 76: BK - BK REV. 77: BK - BK REV. 78: BK - BK REV. 79: BK - BK REV. 80: BK - BK REV. 81: BK - BK REV. 82: BK - BK REV. 83: BK - BK REV. 84: BK - BK REV. 85: BK - BK REV. 86: BK - BK REV. 87: BK - BK REV. 88: BK - BK REV. 89: BK - BK REV. 90: BK - BK REV. 91: BK - BK REV. 92: BK - BK REV. 93: BK - BK REV. 94: BK - BK REV. 95: BK - BK REV. 96: BK - BK REV. 97: BK - BK REV. 98: BK - BK REV. 99: BK - BK REV. 100: BK - BK	REV. 1: BK - BK REV. 2: BK - BK REV. 3: BK - BK REV. 4: BK - BK REV. 5: BK - BK REV. 6: BK - BK REV. 7: BK - BK REV. 8: BK - BK REV. 9: BK - BK REV. 10: BK - BK REV. 11: BK - BK REV. 12: BK - BK REV. 13: BK - BK REV. 14: BK - BK REV. 15: BK - BK REV. 16: BK - BK REV. 17: BK - BK REV. 18: BK - BK REV. 19: BK - BK REV. 20: BK - BK REV. 21: BK - BK REV. 22: BK - BK REV. 23: BK - BK REV. 24: BK - BK REV. 25: BK - BK REV. 26: BK - BK REV. 27: BK - BK REV. 28: BK - BK REV. 29: BK - BK REV. 30: BK - BK REV. 31: BK - BK REV. 32: BK - BK REV. 33: BK - BK REV. 34: BK - BK REV. 35: BK - BK REV. 36: BK - BK REV. 37: BK - BK REV. 38: BK - BK REV. 39: BK - BK REV. 40: BK - BK REV. 41: BK - BK REV. 42: BK - BK REV. 43: BK - BK REV. 44: BK - BK REV. 45: BK - BK REV. 46: BK - BK REV. 47: BK - BK REV. 48: BK - BK REV. 49: BK - BK REV. 50: BK - BK REV. 51: BK - BK REV. 52: BK - BK REV. 53: BK - BK REV. 54: BK - BK REV. 55: BK - BK REV. 56: BK - BK REV. 57: BK - BK REV. 58: BK - BK REV. 59: BK - BK REV. 60: BK - BK REV. 61: BK - BK REV. 62: BK - BK REV. 63: BK - BK REV. 64: BK - BK REV. 65: BK - BK REV. 66: BK - BK REV. 67: BK - BK REV. 68: BK - BK REV. 69: BK - BK REV. 70: BK - BK REV. 71: BK - BK REV. 72: BK - BK REV. 73: BK - BK REV. 74: BK - BK REV. 75: BK - BK REV. 76: BK - BK REV. 77: BK - BK REV. 78: BK - BK REV. 79: BK - BK REV. 80: BK - BK REV. 81: BK - BK REV. 82: BK - BK REV. 83: BK - BK REV. 84: BK - BK REV. 85: BK - BK REV. 86: BK - BK REV. 87: BK - BK REV. 88: BK - BK REV. 89: BK - BK REV. 90: BK - BK REV. 91: BK - BK REV. 92: BK - BK REV. 93: BK - BK REV. 94: BK - BK REV. 95: BK - BK REV. 96: BK - BK REV. 97: BK - BK REV. 98: BK - BK REV. 99: BK - BK REV. 100: BK - BK	<p>PRONTO PRESS-10</p> <p>MAINS 100V - 240V 50/60Hz</p> <p>DATE: 09/10/95 / DFM</p> <p>SCALE: 1:1</p> <p>DATE: Tuesday, September 11, 2001</p>
--	--	--

1. Hydraulic Powerpack
2. Double-acting Hydraulic Cylinder
3.  $\varnothing 5$  Polyamid Tube
4. Pressure Transmitter
5. Lower Valve Manifold
- 5.1 Upper Valve Manifold
6. Relief Valve
7. Accumulator
8. Throttle Plug

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Herewith declares that

<i>Product Name:</i>	ProntoPress-10
<i>Type No:</i>	506
<i>Machine Type:</i>	Mounting Apparatus

is in conformity with the provisions of the following directives:

**Safety of Machinery** 98/37/EEC according to the following standard(s):  
EN ISO 12100-1:2003, EN ISO 12100-2:2003, EN60204-1:1997.


**EMC-Directive** 89/336/EEC and 92/31/EEC according to the following standard(s):  
EN50081-1:1992, EN50082-1:1997.

**Low Voltage Directive** 73/23/EEC according to the following standard(s):  
EN60204-1:1997.

**Supplementary Information** The equipment complies with the American standards:  
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

Date: 01.07.2006

  
Christian Skjold Heyde,  
Vice President, R&D and Production, Struers A/S

**Hersteller** Struers A/S  
Pederstrupvej 84  
DK-2750 Ballerup, Danmark  
Telefon +45 44 600 800

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Date: 01.07.2006

Christian Skjold Heyde,  
Vice- President, R& D et Production, Struers A/S

# Efficient Hot Mounting the Struers Way

## ProntoPress-20

Advanced, programmable mounting press with two cylinders.

Complete process control and built-in database. Synchronous or independent operation on the two mounting cylinders, for highest throughput

## ProntoPress-10

Advanced, programmable single cylinder mounting press.

With complete process control and built-in database

## The Energy Savers

With the newly designed mounting units, half the power produces twice the heat.

Innovative design has merged the heating/cooling coils and the mounting cylinder into one compact and lightweight unit. This allows very rapid heat transfer and therefore very fast heating and cooling. The power consumption of the heating unit has been reduced to 750 W whilst the efficiency, compared to conventional presses, has been increased by about 100%. This means that less than half the energy is required to heat the mounting unit to the set temperature in a shorter time than a conventional press.

## Efficient cooling

Because of the lower mass of the mounting unit, the cooling time is much shorter than that of conventional presses.




Additionally, an "economy setting" during the cooling process reduces the consumption of cooling water to about 20 % of the quantity normally used.

## Recirculation cooling unit

Both ProntoPresses are equipped with a connection for a recirculation cooling unit. That guarantees additional water saving.

## Intermediate ram

By using an intermediate ram, the capacity of the press can be dou-

bled as two mounts can be produced at the same time.

### Easy exchange-six mount sizes

The complete mounting unit can be exchanged in seconds, allowing for easy and fast changes in mount sizes. Six different cylinder diameters are available.

### Perfect Mounting with APC Struers Advanced Process Control takes care of all necessary adjustments

#### Pre-heating

For thermoplastic resins, preheating without or at low pressure will produce better results than otherwise obtainable. On the new Struers presses, both preheating time, temperature and pressure can be adjusted to accommodate specific requirements.

#### Controlled process

For maximum reproducibility the timer can be set to start when the

pre-set temperature is reached. This makes the complete mounting cycle more controlled and reproducible as it eliminates the influence of the initial mounting unit temperature.

#### Pressure or Force

With the ProntoPresses, it is possible for the first time to set and adjust the pressure, instead of the force, exerted on the mount itself. If the cylinder diameter is changed, the press will automatically recalculate the force necessary to obtain the correct pressure. Thus producing every mount similarly, without the need for intermediate readjustments.

#### Controlled cooling, saving of water

Tests have shown that conventional cooling, where water is supplied at a maximum rate, uses far more than necessary. Also, different resins require different cooling gradients, i.e. thermoplastic resins need far longer cooling times than thermosetting ones.

With ProntoPress, however, three different rates of cooling are available, enabling the cooling rate to be more accurately set for different resins. Instead of reducing the amount of water by manually adjusting the flow at the source, the ProntoPresses have preprogrammed levels delivering pulses of cooling water at set intervals.

This results in almost the same cooling rate as a constant flow of water, whilst at the same time reducing the amount of water used by up to 80%.

#### Setting your own user interface


Due to the advanced software it is very easy to change the screen prompts, choosing between four of the major languages.

Also physical units can be displayed according to your normal usage. Do you prefer to display the parameters as pressure (psi, bar, MPa) or force (Lbf, kp, kN)?

Are you used to expressing tem-

*Prontopress-20 with its two cylinders is indispensable when mounting large series of samples in different diameters*

perature in Celsius or Fahrenheit? These settings are easily configured when using the machine for the first time.

## How to remember all mounting data

### Database incorporated

To store all mounting data for different materials or resins, a database of parameters has been incorporated. ProntoPress-10 can save up to 100 mounting methods, whereas ProntoPress-20 can save up to 200 mounting methods, 100 per cylinder. Struers methods, included in both presses, contain 10 Struers methods with the default values for the resins and relevant mounting conditions. The press will automatically set the correct parameters for the actual cylinder diameter.

For mounts where transparency is a must, Specifast is used

The use of an electro-hydraulic pump guarantees smooth movement of the ram. Even very thin samples can be placed and will remain in position when the ram is moved up or down.

## New top closure

Also a new top closure has been developed, combining the advantages of both bayonet and threaded closure. Equipped with a triple external thread, the new closure will always “catch” like a bayonet and the threads will force the upper ram down into the cylinder, highly facilitating the work for the operator. As the threads are external, contamination with resin is not possible.

## Elevator with steering arm and built in safety switch

The top closure is positioned in an elevator with steering arm. When a mounting process is finished, the top closure is unscrewed. The elevator then raises the top closure, allowing it to be swung to the side, rendering the mounting cylinder easily accessible.

To mount the next sample, simply turn the elevator back into the correct position and press the top closure down into the cylinder. If the elevator is not pressed down, a safety switch will not allow the press to start.

On the display both the information about the methods location in the database, the resin and the mounting parameters are shown

## Technical Data

	ProntoPress-20 PRO20	ProntoPress-10 PRO10
<b>Code</b>		
<b>Mounting Unit</b>		
No. of units	2	1
Mounting dia. interchangeable	25, 30, 40, 50 mm and 1 1/4" and 1 1/2"	
Dual mounts	Optional intermediate ram	
<b>Mounting Parameters</b>		
Force range, kN adjustable in steps of	0 - 50 5	
Pressure range, MPa adjustable in steps of (depending on the diameter)	0 - 25 to 100 2.5 to 10	
Temperature range, degr. C adjustable in steps of	100 - 180 10	
Time for preheating, heating and cooling, min adjustable in steps of	0 - 15:00 0:30	
Cooling rate	Economy: Low: High:	20% of full flow 2% of full flow in the first 3 min. then 20% of full flow Full flow
<b>Software and Electronics</b>		
LC display	4 x 40 characters	
Controls	Touch pads	
Memory	EPROM, RAM and non-volatile RAM 8 kb	
Database	2 x 100 mounting methods Copying and numeric/alphanumeric editing facilities	100 mounting methods
Process control for pre-heating and heating	Time from start / time after temperature is reached Temperature Force/pressure	
Process control for cooling	Time Cooling rate	
Manual control	Heating: Cooling	Temperature Time, cooling rate
Read out during mounting process	Actual process step for both cylinders Remaining process time Preset/actual force/pressure Cooling rate	
<b>Safety Standards</b> (conforming to: CE and FCC)		
General	89/392EEC & 91/386/EEC & 93/44/EEC EN292-1, EN292-2, EN60204-1 (IEC 204-1)	
EMC	89/336/EEC & 92/231/EEC EN50081-2, EN50082-2, FCC part 15, subpart B	
<b>Supply</b>		
Electrical, no. of phases	1	
Power consumption, W		
Max. total	1500	1500
Mounting unit	2 x 750	750
Water	Tap water: min. 1 bar, max. 10 bar	
Recirculating cooling unit	Automatic start/stop of TRECI connected with ROPRE	
<b>Dimensions and Weight</b>		
Width x depth x height, mm (Mounting unit with top closure installed)	420 x 580 x 420	
Weight, kg	51	38

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup, Denmark  
Phone +45 44 600 800  
Fax +45 44 600 801  
struers@struers.dk  
www.struers.com

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**Struers Ltd.**
**Struers Inc.**