RotoForce-1

Instruction Manual

Manual No.: 15227001

Date of Release 20.02.1998

RotoForce-1
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide	15
Quick Reference Guide	23

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 1998.

Struers A/S
Valhøjs Allé 176
DK-2610 Rødovre/Copenhagen
Denmark
Telephone +45 36 70 35 00
Fax +45 38 27 27 01

RotoForce-1

Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the machine according to the Instruction Manual. The machine must be placed in an adequate working position.
2. Be sure that the actual voltage corresponds to the voltage stated on the back of the machine. The machine must be earthed.
3. Keep your hands clear of the supporting plate or specimen mover plate when lowering RotoForce-1.
4. Do not touch the rotating parts during operation.
5. If you observe malfunctions or hear unusual noises - stop the machine and call technical service.
6. Blue lubricant: follow the current safety rules for handling, mixing, filling, emptying and disposal of the alcohol-based lubricant.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of Packing	2
Getting Acquainted with RotoForce-1	3
Setting Up RotoForce-1 on RotoPol-11	4
RotoForce-1 on RotoPol-21/22/25	4
Disc Margin	5
Connecting RotoForce-1	5
Supplying Power	5
Electrical Connection	5
Connection to Compressed Air	5
Connection to the Preparation Equipment	6
RotoPol alone	6
RotoPol with Multidoser and/or RotoCom	6
Choosing between Menus	6
Setting the Language	7
Mounting Lupo	8
2. Operation	
Using the Controls	9
Front Panel Controls of RotoForce-1 and Lupo	9
Front Panel Controls	10
RotoForce-1	10
Lupo	10
The Handle	11
Changing to Another Lupo Unit (Accessory)	11
Changing Lubricant Type / Refilling	11
Running a Preparation Method	12
Using a Specimen Mover plate	12
Using the Specimen Holder	12
Adjusting the Force	12
Adjusting the Preparation Time	12
Operation of Lupo	13
Operation of RotoForce-1	13

1. Getting Started

Checking the Contents of Packing

In the packing box you should find the following parts:


- 1 RotoForce-1
- 1 Mains cable
- 1 Lupo automatic drip lubricator
- 1 Nozzle arm (to be used with Multidoser)
- 1 Pressure hose, 1/4" 3 m
- 1 Nylon tube ø8x6
- 1 Hose clamp
- 1 Quick coupling
- 1 Hose nipple
- 6 Screws M4x10
- 2 Screw M5x12
- 1 Screw M3x3
- 3 Washers ø16.5x30x3
- 3 Washers ø16.5x30x1
- 1 Nut M16
- 1 Socket spanner
- 1 Hexagon key 1.5
- 1 Hexagon key 2.5
- 3 Instruction manuals
- 1 Spare Part and Diagrams manual

Remove RotoForce-1 from the transportation box.

RotoForce-1
Instruction Manual

**Getting Acquainted with
RotoForce-1**

Take a moment to familiarise yourself with the location and names of the RotoForce-1 components.

RotoForce-1 on RotoPol-11 (optional)

- A** Control keys and display.
- B** Lupo drip lubricator.
- C** Release handle.
- D** Support column.
- E** Pressure feet.
- F** Fastening screws and support hole.

**Setting Up RotoForce-1
on RotoPol-11**


- Remove the cap covering the support hole placed on the top of the RotoPol-11 cabinet.
- Remove the two caps covering the fastening screws on the left-hand side of the RotoPol cabinet by lifting them off with a knife.
- Lead the supporting column of RotoForce-1 down into the support hole in RotoPol-11. Turn the column until the specimen mover plate is positioned correctly over the preparation disc. (see figure)
- Mount a preparation disc on RotoPol-11 and adjust the distance between disc and specimen mover plate to about 1 - 1.5 mm.
- Tighten the two fastening screws on the left-hand side of the RotoPol-11 cabinet and mount the two caps again.

**RotoForce-1
on RotoPol-21/22/25**

Place two working tables of equal height at a distance about half the length of RotoPol. Balance RotoPol between the tables.

- Lead the supporting column of RotoForce-1 down into the support hole of RotoPol.
- Turn the column until the specimen mover plate is positioned correctly over the preparation disc.
- Fasten the nut and washer from the transportation box loosely under the column of RotoForce-1 below RotoPol.
- Adjust the height of the specimen mover in relation to the preparation disc:
 - Mount a preparation disc on RotoPol.
 - Mount a specimen mover plate on RotoForce-1.
 - Lower RotoForce-1 into working position.
 - Measure the distance between the specimen mover plate and the preparation disc.
 - The correct distance should be 1 - 1.5 mm. If the distance is not correct, remove RotoForce-1 again.
 - Grease 1 or more extra washers (1-3 mm) and place them on the thread at the end of the column. (the grease makes the washers stick to the column.)
 - Put RotoForce-1 back onto RotoPol and make sure the distance is correct now.
 - Fasten the bottom nut tightly.

Disc Margin

Before final tightening of the two fastening screws or the bottom nut you must check the margin between the specimen holder plate hole and the edge of the grinding / polishing disc. This distance should be:

Diameter of the Holes	Disc Margin
40 mm	1.5 mm
1½"	2.5 mm
1¼"	6.0 mm
30 mm	6.5 mm
25 mm	9.0 mm

IMPORTANT

Having adjusted the disc margin, tighten the two fastening screws or the bottom nut firmly. The column should now remain in position.

**Connecting RotoForce-1
Supplying Power**

RotoForce-1 is supplied with electric cable, air hose and connection cable.

IMPORTANT

Check that the mains voltage corresponds to the voltage stated on the type plate on the back of the machine.

Electrical Connection

The mains cable of RotoPol is either equipped with a standard electric male plug, type "F" or with a male flat blade attachment plug. If the type of plug delivered with the equipment has not been approved for local use, the plug must be replaced with an approved type. The leads have to be connected as following:

- Yellow/green: earth
 - Brown: phase
 - Blue or black: neutral
- Insert the female plug of the power cable into the socket on the back of the machine.
 - Connect the male plug to a grounded electric outlet.
- Mount an air hose on the enclosed quick coupling and secure with the enclosed hose clamp.
 - Connect the quick coupling to the factory-mounted air hose of RotoForce-1.

Connection to Compressed Air

IMPORTANT

The air pressure must be between 5 bar (73 psi) and 6.5 bar (95 psi)

RotoForce-1
Instruction Manual

**Connection to
the Preparation Equipment**
RotoPol alone


Connect the communication cable to the back of RotoPol (see the Instruction Manual of RotoPol).

*RotoPol with Multidoser and/or
RotoCom*

Connect the communication cable to the back of Multidoser (see the Instruction Manual of Multidoser)

Choosing between Menus

Press the MENU  key to move from one screen to another:

Setting the Language

30N 01:00
STOP

LANGUAGE
ENGLISH

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
30N 01:00
STOP

Press MENU  to activate the set-up menu.

The LANGUAGE screen appears.

Press SELECT  to choose English, German, French or Japanese.

Press STOP  to accept the choice. The Operation Menu will appear again in the language you have chosen.

Remember to set the same language on all the connected RotoSystem modules.


Mounting Lupo

- Remove Lupo from its box
- Guide Lupo onto the connector on the supporting arm on RotoForce-1.
- Remove the top lid of the bottle and fill with lubricant.
- Remount the top lid, and Lupo is ready for use.


2. Operation

Using the Controls
Front Panel Controls of
RotoForce-1 and Lupo

RotoForce-1





Lupo



RotoForce-1
Instruction Manual

Front Panel Controls

RotoForce-1

Name	Key	Function	Name	Key	Function
① MAIN SWITCH		The main switch is located on the back of RotoForce-1.	⑤ MENU		Activates the set-up menu or moves down in the menu structure.
② AUTO		Turns Auto mode on or off.	⑥ SELECT		Changes values in the set-up menu.
③ START		Starts RotoForce-1.	⑦ ⑧ FORCE		Select FORCE: ▲ Increases the force ▼ Decreases the force
④ STOP		Stops RotoForce-1. Exits the set-up menu.	⑨ ⑩ TIME		Select TIME: ▲ Increases the time ▼ Decreases the time

Lupo

Name	Key	Function	Name	Key	Function
⑪ DOSING SWITCH		Switches between different step dosing levels.	START		Manual predosing.
⑫ CONTROL SWITCH		Switches between Auto, Start and Stop mode. Lupo is supplied with power direct from RotoForce-1 and does not need a separate power switch.	STOP		Stop the predosing.
			AUTO	Auto	When Auto mode is active start and stop of Lupo can be controlled from RotoForce-1. When RotoForce-1 is started/stopped Lupo also starts/stops.

The Handle

The handle of RotoForce-1 is used to move RotoForce-1 into an upright position. In this position it is possible to exchange specimen mover plates and preparation discs. As RotoForce-1 is equipped with a gas spring, it will be moved upwards automatically after the handle is moved backwards.

To move RotoForce-1 down again, simply press it down until the locking pins engage in the bottom plate.

Changing to Another Lupo Unit (Accessory)

- Set the mode switch on Lupo to STOP  position.
- Remove the current Lupo.
- Place the other Lupo (accessory) with an alternative lubricant in the connector on the supporting arm.

Changing Lubricant Type / Refilling

- Set the mode switch on Lupo to STOP  position.
 - Remove Lupo from RotoForce-1.
 - Hold the Lupo lubricant bottle firmly and remove the top lid.
 - Empty the bottle, still mounted on Lupo, of any remaining lubricant.
 - Fill the bottle with mild, soapy water.
 - Put Lupo back on the supporting arm.
 - Set the mode switch on Lupo to START .
 - Set the dosing control to max. dosing (step 10) and let the water run through the nozzle for 2-3-min.
- Repeat the above points but fill with tap water to rinse out any soap residue.
- Refill with lubricant and remount the top lid.

IMPORTANT

Always avoid turning the socket cap while unscrewing the Lupo bottle, as the turning movement may cause damage to the socket cap gasket.

Running a Preparation Method

Using a Specimen Mover Plate


- Move RotoForce-1 into the upright position.
- Remove the three screws securing the mover plate or slide the specimen holder off the legs.
- Put a specimen mover plate in position and secure it with the three screws.

Using the Specimen Holder

- Move RotoForce-1 into the upright position.
- Remove the three screws securing the specimen mover plate.
- Slide the specimen holder evenly onto the three legs.
- Insert the correct preparation disc and press RotoForce-1 down into the correct position.

Adjusting the Force


Press FORCE (F)▲ or (F)▼ to adjust the force to the desired value.

The force displayed is the force applied per specimen. Also when working with the specimen holder ROFLU, as the three pressure feet will transfer the force to the specimen holder where three specimens are clamped. It is not necessary to apply higher forces.

Adjusting the Preparation Time

Press TIME (T)▲ or (T)▼ to adjust the preparation time.

Operation of Lupo

- Adjust the dosing level
- Set Lupo to Auto mode

Operation of RotoForce-1

- After setting the correct time and force on RotoForce-1 and the dosing level on Lupo, press Auto on RotoForce-1.
- Press START ◊ on the RotoPol grinding/polishing machine to start both machines simultaneously.
- When the pre-set time has elapsed, all machines will stop simultaneously.

RotoForce-1
Instruction Manual

Reference Guide

Table of Contents	Page
1. Force Reduction	16
2 Grinding	
Plane Grinding.....	17
Fine Grinding.....	17
Recommendations.....	17
3. Accessories	18
4. Trouble-shooting	19
5. Maintenance	
Daily.....	20
Weekly.....	20
Lupo.....	20
6. Technical Data	21


1. Force Reduction


RotoForce-1 is equipped for automatic reduction of the force towards the end of the preparation.



If FORCE REDUCTION is enabled, the force will automatically be reduced to 50% of the pre-set value for the last third of the preparation step. This will result in lower material removal and smaller scratches than usual. The time for the subsequent preparation steps can be reduced as the amount of damaged material to be removed will be reduced.


Press MENU  until you reach the FORCE REDUCTION screen.

Press SELECT  to toggle between NO and YES.

Press STOP  to accept the choice. Operation Menu will appear again.

2 Grinding

The grinding process can be divided into plane grinding **PG** and fine grinding **FG**. In the following short descriptions you will find a number of useful tips for grinding with RotoForce-1.

Plane Grinding

PG is unnecessary when the specimens are prepared as single specimens, as all such specimens are levelled individually. After sectioning, using a Struers wet-cutting machine, the surface of the specimens is plane enough to proceed directly to a fine grinding stage.

When using the specimen holder ROFLU the specimens must be ground to the same level using a **PG** step.

Fine Grinding

FG is used to remove the deformation introduced during **PG** or abrasive wet-cutting, and to produce specimens with a surface ready for polishing. We recommend one-step **FG** like MD-Allegro, MD/DP-Plan or MD/DP-Pan.

Recommendations

Do not use RotoForce-1 with specimen mover plates for plane grinding. **PG** is normally unnecessary, and the use of coarse abrasives can result in unplane specimens.

If, for whatever reason, it is necessary to grind using coarse abrasive when working with specimen mover plates, the planeness may be improved by carrying out the following:

- Check that the distance between specimen mover plate and preparation disc is maximum 1 - 1.5 mm.
- The height of the specimen should not exceed 0.7 x specimen diameter. Example: a specimen with a diameter of 30 mm should not be higher than $30 \times 0.7 = 21$ mm.
- Use a grain size as small as possible.
- Use a mounting resin with a wear resistance similar to the specimens' wear resistance.
- Use 150 rpm speed for both the grinding disc and specimen mover.
- Use low force.


3. Accessories

Specification	Code
<i>Specimen mover plates</i> for 3 specimens, Ø25 mm/1" for 3 specimens, Ø30 mm for 3 specimens, Ø13" for 3 specimens, Ø12" for 3 specimens, Ø40 mm Without holes	ROFIN ROFTH ROFQU ROFHA ROFYF ROFNO
<i>Specimen holder</i> for clamping of 3 specimens 15 - 40 mm dia.	ROFLU
<i>Levelling device</i> for levelling of specimens in specimen holders	ROFFU
LUPO, optional extra drip lubricator	PEDAU

Remember...

Struers offers a comprehensive range of consumables for both grinding and polishing. Please ask for separate leaflets.

4. Trouble-shooting

Display/Error	Cause	Action
RotoForce-1		
OVERLOAD	RotoForce-1 is stopped due to overloading of the motor.	Wait for the motor to cool down and the PAUSE/OVERLOAD message will appear on the display.
PAUSE/ OVERLOAD	The motor has cooled down after overloading and RotoForce-1 is ready for use.	Restart RotoForce-1.
PAUSE/ NO AIR	Air pressure too low.	Supply correct air pressure and restart
Lupo		
Bad or missing dosing	Dosing nozzle clogged	Clean the nozzle with alcohol
Lubricant leaking below the bottle	The bottle has not been fastened correctly	Fasten tightly. Remember to hold back the socket cap itself to avoid damage to the bottom gasket
Preparation Error		
Unplane specimens	Abrasive too coarse	See section Grinding.
Japanese Texts	RotoForce-1 has the possibility of showing the screen texts in Japanese. When choosing a European language in the LANGUAGE menu, you might accidentally press  MENU at the Japanese option.	To return to English, German or French texts, just follow this procedure: <ul style="list-style-type: none"> ■ Press STOP . ■ Press MENU  two times. ■ Press SELECT  until you reach the language you prefer. ■ Press STOP .

5. Maintenance

Daily

Clean all accessible surfaces of the machine with a moist cloth

Weekly

Lupo

To avoid clogging of the Lupo nozzle, you should clean the nozzle once a week:

- Set the mode switch on Lupo to STOP  position.
- Remove Lupo from the supporting arm.
- Hold the Lupo lubricant bottle firmly and remove the top lid.
- Empty the bottle, still mounted on Lupo, of any remaining lubricant.
 - Fill the bottle with mild, soapy water.
 - Put Lupo back on the supporting arm.
 - Set the mode switch on Lupo to START .
 - Set the dosing control to max. dosing (step 10) and let the water run through the nozzle for 2-3-min.
 - Repeat the above points but fill with tap water to rinse out any soap residue.
 - Refill with lubricant and remount the top lid.
 - Clean the O-rings and lubricate them using High Vacuum Grease

6. Technical Data

Subject		Specifications		
		Metric/International	US	
Disc	Rotational speed	150 rpm	150 rpm	
	Rotational direction	Counter clockwise	Counter clockwise	
	Motor output 50/60Hz	45W	45W	
	Torque at disc	Continuous 50/60HZ	1.25Nm	1.25Nm
		Max. 50/60HZ	1.7Nm	1.7Nm
	Air pressure	Min. 5.0 max. 6.5 bar	Min. 73 max. 95 psi	
Specimens	Force on single specimens	5-50N in 5N steps	5-50N in 5N steps	
	Specimens in specimen mover plate	Ø25-40 mm	1-1 1/2"	
	Specimen height	8-28 mm	0.3-1.1"	
Dimensions and Weight	Width	200 mm	7.9"	
	Depth	350 mm	13.8"	
	Height	445 mm	17.5"	
	Weight	11.2 kg	24.7 lb	

RotoForce-1
Instruction Manual

Subject		Specifications
Software and electronics	LC display	2x16 characters
	Controls	Touch pad
	Memory	EPROM, RAM and non-volatile RAM
	Automatic start/stop	With: Lupo
	Communication including automatic start/stop	With: RotoPol-11/21/22/25 Multidoser
	Data exchange	With: RotoCom
Environment	Safety Standards	General 89/392/EEC & 91/386/EEC & 93/44/EEC EN292-1, EN292-2, EN60204-1 (IEC 204-1) EMC 89/336/EEC & 92/231/EEC EN50081-1, EN50082-1 FCC part 15, subpart B
	Noise Level	Approx. 52 dB (A) at idle running, at a distance of 1.0 m/39.4" from the machine
	Surrounding temperature	5-40°C/41-104°F
	Humidity	0-95% RH non condensing
Supply Voltage	Power consumption continued	90 W
	No. of phases	1
	Voltage / frequency:	Max. continuous load
	100V/50Hz	0.6A
	100-120V/60Hz	0.5A
220-240V/50Hz	0.2A	
220-240V/60Hz	0.22A	

Quick Reference Card

Running a Preparation Method

Using a Specimen Mover plate

- Move RotoForce-1 into the upright position.
 - Remove the three screws securing the mover plate or slide the specimen holder of the legs.
 - Put a specimen mover plate in position and secure it with the three screws.

Using the Specimen Holder

- Move RotoForce-1 into the upright position.
- Remove the three screws securing the specimen mover plate.
- Slide the specimen holder evenly onto the three legs.
- Insert the correct preparation disc and press RotoForce-1 down into the correct position.

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
Struers A/S
Valhøjs Allé 176
DK-2610 Rødovre, Denmark
Telephone 36 70 35 00

Herewith declares that

<i>Product Name:</i>	RotoForce-1
<i>Type No:</i>	522
<i>Machine Type:</i>	Specimen mover

is in conformity with the provisions of the following directives:


Safety of Machinery 89/392/EEC and 91/368/EEC and 93/44/EEC according to the following standard(s):
EN292-1:1991, EN292-2:1991, EN60204-1:1992.

EMC-Directive 89/336/EEC and 92/31/EEC according to the following standard(s):
EN50081-1:1991, EN50082-1:1992.

Low Voltage Directive 73/23/EEC and 93/68/EEC according to the following standard(s):
EN60204-1:1992.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

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Klaus Kisbøll, R&D Manager, Struers A/S
Hersteller Struers A/S
Klaus Kisbøll, Entwicklungsleiter, Struers A/S
Fabricant Struers A/S
Klaus Kisbøll, Chef du service developpement, Struers A/S